

PERKIN ELMER CETUS

REFERENCES

1. Saiki, R.K., Gelfand, D.H., Stoffel, S., Scharf, S.J., Higuchi, R., Horn, G.T., Mullis, K.B., and Erlich, H.A. 1986. Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase. *Science* 239: 487-491.
2. Saiki, R.K., Scharf, S., Faloona, F., Mullis, K.B., Horn, G.T., Erlich, H.A., and Arnheim, N. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* 230: 1350-1354.
3. Scharf, S.J., Horn, G.T., and Erlich, H.A. 1986. Direct cloning and sequence analysis of enzymatically amplified genomic sequences. *Science* 233: 1076-1078.
4. Gyllenstein, U.B. and Erlich, H.A. 1988. Generation of single stranded DNA by the polymerase chain reaction and its application to direct sequencing of the HLA DQ alpha locus. *Proc. Natl. Acad. Sci. USA*. 85: 7652-7656.
5. Kawasaki, E.S., Clark, S.C., Coyne, M.Y., Smith, S.D., Champlin, R., White, O.N., and McCormick, F.P. 1988. Diagnosis of chronic myeloid and acute lymphocytic leukemias by detection of leukemia-specific mRNA sequences amplified *in vitro*. *Proc. Natl. Acad. Sci. USA*. 85: 5698-5702.
6. Lawyer, F.C., Stoffel, S., Saiki, R.K., Myambo, K.B., Drummond, R., and Gelfand, D.H. 1989. Isolation, characterization, and expression in *E. coli* of the DNA Polymerase gene from the extreme thermophile, *Thermus aquaticus*. *J. Biol. Chem.* 264: 6427-6437.
7. Innis, M.A., Myambo, K.B., Gelfand, D.H., and Brow, M.A.D. 1988. DNA sequencing with *Thermus aquaticus* DNA polymerase. *Proc. Natl. Acad. Sci. USA*. 85: 9436-9440.
8. Higuchi, R. and Kwok, S. 1989. Avoiding false positives with PCR. *Nature* 339: 237-238.
9. Gelfand, D.H. Cetus Corp. 1988. Personal communication.
10. McConlogue, L., Brow, M.A.D., and Innis, M.A. 1988. Structure-independent DNA amplification by PCR using 7-deaza-2'-deoxyguanosine. *Nucleic Acids Research* 16: 9869.
11. Livingston, D.M., Hinkle, D.C., and Richardson, C.C. 1975. Deoxyribonucleic acid polymerase III of *Escherichia coli*: purification and properties. *J. Biol. Chem.* 250: 461-469.
12. Sanger, F., Coulson, A.R., Hong, G.F., Hill, D.F., and Petersen, G.B. 1982. Nucleotide sequence of bacteriophage lambda DNA. *J. Mol. Biol.* 162: 729.
13. Mullis, K.B. and Faloona, F.A. 1987. Specific synthesis of DNA *in vitro* via a polymerase-catalyzed chain reaction. In: Wu, R., ed. *Methods in Enzymology*, Vol. 155. San Diego, CA, Academic Press pp. 335-350.
14. Jeffreys, A.J., Wilson, V., Neumann, R., and Keyte, J., 1988. Amplification of human minisatellites by the polymerase chain reaction: towards DNA fingerprinting of single cells. *Nucleic Acids Research*, 16: 10953-10971.
15. Gelfand, D.H., and White, T.S. 1990. Thermostable DNA Polymerases. In: Innis M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. *PCR Protocols, A Guide to Methods and Applications*, Academic Press, Inc., San Diego, CA, pp 129-141.
16. Innis M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. 1990. *PCR Protocols, A Guide to Methods and Applications*. Academic Press, Inc., San Diego, CA.
17. Kwok, S. 1990. Chapter 17. Procedures to minimize PCR-product carry-over, and Orrego, C. 1990. Chapter 54. Organizing a laboratory for PCR work. In: Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. 1990. *PCR Protocols, A Guide to Methods and Applications*. Academic Press, Inc., San Diego, CA.
18. Rychlik, W., Spencer, W.J., and Rhoads, R.E. 1990. Optimization of the annealing temperature for DNA amplification *in vitro*. *Nuc. Acids Res.* 18: 6409-6412.
19. Abramson, R.D., et al. Cetus Corporation (Manuscript in Preparation).

PERKIN ELMER CETUS

ORDERING INFORMATION

**PRODUCTS ARE NOT FOR MEDICAL OR DIAGNOSTIC USE
FOR RESEARCH USE ONLY**

<p><i>In the United States:</i></p> <p>For PCR Technical support call toll-free 1-800-762-4001. FAX: 203-834-4391</p> <p>For easy ordering from PE XPRESS call toll-free 1-800-762-4002.</p> <p>To inquire about status of orders or to request a catalog from PE XPRESS call toll-free 1-800-762-4003.</p> <p>For product literature on the full line of Perkin-Elmer Cetus Instruments and Reagents call toll-free 1-800-762-4000.</p>	<p><i>Outside of the United States:</i></p> <p>For the number of your local Perkin-Elmer sales representative call:</p> <table> <tr><td>Austria</td><td>(222) 623 1010</td></tr> <tr><td>Belgium</td><td>02/21 20 734</td></tr> <tr><td>Canada</td><td>514-737-7575</td></tr> <tr><td>Denmark</td><td>42/811555</td></tr> <tr><td>Far East-Australia</td><td>3-560-4566</td></tr> <tr><td>Finland</td><td>90/880 144</td></tr> <tr><td>France</td><td>(1) 30 85 63 63</td></tr> <tr><td>Germany</td><td>06103-708-118</td></tr> <tr><td>Italy</td><td>039/83841</td></tr> <tr><td>Latin America-Mexico</td><td>5-651-7077</td></tr> <tr><td>Netherlands</td><td>0 1820-75481</td></tr> <tr><td>Norway</td><td>02/686565</td></tr> <tr><td>Portugal</td><td>1762172</td></tr> <tr><td>Spain</td><td>94/447-1021</td></tr> <tr><td>Sweden</td><td>08/7330010</td></tr> <tr><td>Switzerland</td><td>061/6925740</td></tr> <tr><td>United Kingdom</td><td>(02407)4411</td></tr> </table>	Austria	(222) 623 1010	Belgium	02/21 20 734	Canada	514-737-7575	Denmark	42/811555	Far East-Australia	3-560-4566	Finland	90/880 144	France	(1) 30 85 63 63	Germany	06103-708-118	Italy	039/83841	Latin America-Mexico	5-651-7077	Netherlands	0 1820-75481	Norway	02/686565	Portugal	1762172	Spain	94/447-1021	Sweden	08/7330010	Switzerland	061/6925740	United Kingdom	(02407)4411
Austria	(222) 623 1010																																		
Belgium	02/21 20 734																																		
Canada	514-737-7575																																		
Denmark	42/811555																																		
Far East-Australia	3-560-4566																																		
Finland	90/880 144																																		
France	(1) 30 85 63 63																																		
Germany	06103-708-118																																		
Italy	039/83841																																		
Latin America-Mexico	5-651-7077																																		
Netherlands	0 1820-75481																																		
Norway	02/686565																																		
Portugal	1762172																																		
Spain	94/447-1021																																		
Sweden	08/7330010																																		
Switzerland	061/6925740																																		
United Kingdom	(02407)4411																																		
<p><i>In Europe:</i></p> <p>For PCR Technical support call Perkin-Elmer's European Sales Support Center in Vaterstetten, Germany at 49/8106/381-112. FAX: 49/8106/6697</p>																																			

GeneAmp™ PCR Instrument Systems from Perkin-Elmer Cetus Instruments:

- N801-0150 DNA Thermal Cycler, 120 V, 60 Hz
- N801-0177 DNA Thermal Cycler, 240 V, 50 Hz
- N801-0189 DNA Thermal Cycler, 100 V, 50/60 Hz

- N801-0100 DNA Thermal Cycler 480, 120 V, 60 Hz
- N801-0101 DNA Thermal Cycler 480, 240 V, 50 Hz
- N801-0102 DNA Thermal Cycler 480, 100 V, 50/60 Hz

- N801-0001 GeneAmp PCR System 9600, 120 V, 60 Hz
- N801-0002 GeneAmp PCR System 9600, 240 V, 50 Hz
- N801-0003 GeneAmp PCR System 9600, 100 V, 50/60 Hz

Perkin-Elmer® is a registered trademark of The Perkin-Elmer Corporation.
Cetus® is a registered trademark of the Cetus Corporation.
GeneAmp®, AmpliTaq®, and GeneAmp™ are registered trademarks/trademark of Cetus Corporation, licensed to Perkin-Elmer Cetus
Instruments.

Tween 20® is a registered trademark of ICI.
NuSieve® and SeaKem® are registered trademarks of FMC Corporation.

Copyright 1991 by Perkin-Elmer Cetus Instruments. All rights reserved.

BIO-23
55635-5/91

PERKIN ELMER CETUS
761 Main Ave., Norwalk, CT 06859

GeneAmp[®] HIV-1 Reagents

**NOT FOR DIAGNOSTIC OR MEDICAL USE
FOR RESEARCH USE ONLY**
See notice to purchaser

Part No. N808-0015

The Perkin-Elmer Cetus GeneAmp[®] HIV-1 Reagents are designed for use in detecting the presence of HIV-1 gag gene sequences in experimental samples using the GeneAmp[®] Polymerase Chain Reaction (PCR) process. These reagents include both positive and negative HIV-1 control DNA templates. The positive control is plasmid DNA that contains the entire genome of the HIV26 isolate. The HIV26 genome has been rearranged to interrupt the *pol* gene region and thus block infectivity. The plasmid has been shown to be non-infectious in cell culture.¹ Human placental DNA that is free of HIV-1 sequences is provided as a negative control. Included are primers SK38 and SK39 which can be used to amplify a 115 bp sequence in the gag region of the HIV-1 genome.¹ Detection of the 115 bp PCR product can be achieved using the SK19 probe after suitable labeling. The GeneAmp[®] HIV-1 Reagents are optimized for use with the GeneAmp[®] PCR Core Reagents (Part No. N808-0009) or the GeneAmp PCR Reagent Kits with Native Taq DNA Polymerase (Part No. N801-0043) or AmpliTaq[®] DNA Polymerase (Part No. N801-0055). Sufficient primer and probe reagents are provided for 100 amplification reactions of 100 μ L each. Depending upon the dilutions used, sufficient control materials are provided for 10-100 amplification reactions.

The GeneAmp PCR process is a simple and powerful method¹, invented by K. Mullis and patented by Cetus Corporation, which allows *in vitro* amplification of DNA segments¹ (including cDNA)¹ through a succession of incubation steps at different temperatures. Typically, the double stranded DNA is heat denatured, the two primers complementary to the 3' boundaries of the target segment are annealed at a low temperature and then extended at an intermediate temperature. One set of these three consecutive steps is referred to as a cycle. For the GeneAmp[®] HIV-1 Reagents it is possible to simultaneously anneal and extend at the intermediate temperature thereby reducing the number of steps by one, leading to a decrease in cycling time while still maintaining specificity and the level of amplification.

LIST OF KIT COMPONENTS

Store the GeneAmp[®] HIV-1 Reagents at -20°C in a constant temperature freezer. If stored under the proper conditions, the product will maintain performance until the control date printed on the label.

Reagent	Amount	Volume	Concentration	Comments
HIV-1 Positive Control DNA	10 ³ copies	100 μ L	10 ³ copies/ μ L	In 10 μ g/mL human placental DNA, 1 mM EDTA, 10 mM NaCl, 10 mM Tris-HCl, pH 8.0.
HIV-1 Negative Control DNA from human placenta	1 μ g	100 μ L	10 μ g/mL	In 1 mM EDTA, 10 mM NaCl, 10 mM Tris-HCl, pH 8.0.
HIV-1 Primer 1 SK38	5 nmoles	200 μ L	25 μ M	In 10 mM Tris-HCl, pH 8.3.
HIV-1 Primer 2 SK39	5 nmoles	200 μ L	25 μ M	In 10 mM Tris-HCl, pH 8.3.
HIV-1 Probe SK19	100 pmoles	50 μ L	2 μ M	In 10 mM Tris-HCl, pH 8.3.

REAGENTS NEEDED BUT SOLD SEPARATELY

Reagent	Part No.
Ampli [®] Taq [®] DNA Polymerase*	N801-0060
GeneAmp dNTPs*	N808-0007
GeneAmp 10X PCR Buffer II & MgCl ₂ Solution*	N808-0010
*These items can be purchased together by ordering the GeneAmp PCR Core Reagents	N808-0009

Notice to Purchaser: Disclaimer of License

This product is designed for use in the GeneAmp[®] Polymerase Chain Reaction (PCR) process covered by patents issued to Cetus Corporation, but no license under these patents to use the GeneAmp PCR process is conveyed by the purchase of this product. A license to use the GeneAmp PCR process for certain research and development activities accompanies the purchase and use of Perkin-Elmer Cetus' GeneAmp[®] PCR Core Reagents, GeneAmp[®] PCR Reagent Kits, and certain other products, or is available from Cetus Corporation. For information contact Cetus Director of Business Development, PCR Division, at 1400 Fifty-Third Street, Emeryville, CA 94608. Tel. (415) 420-3300.

ATTORNEY'S EYES ONLY
RESTRICTED

RMS 71151

PERKIN ELMER CETUS

PROTOCOL FOR DNA AMPLIFICATION

I. General Advice

Due to the enormous amplification possible with the GeneAmp[®] PCR process, small levels of DNA contamination, especially from previous PCR amplification reactions and positive control templates, can result in product even in the absence of purposefully added template DNA.* If possible, all reactions should be set up in an area separate from PCR product analysis or sample preparation. Dedicated pipettes, vessels and solutions for DNA preparation, reaction mixing, and sample analysis will minimize cross contamination.

A master mix of reagents (water, buffer, dNTPs, primers and enzyme) for all samples can be prepared first, then aliquoted to individual tubes. Magnesium chloride and the template DNA are then added. Using such mixes will minimize reagent pipetting losses, increase accuracy, and reduce the number of reagent transfers.

Amplifications are performed in 100 μ L of reaction mix in capped 0.5 mL polypropylene microcentrifuge tubes. Perkin-Elmer Cetus 0.5 mL GeneAmp[®] Reaction Tubes (Part No. N801-0180) provide the best heat transfer when using the Perkin-Elmer Cetus DNA Thermal Cycler, because of their uniform fit in the wells. Since DNA may stick to plastic and since nucleases are often found on surfaces, it may be preferable to use sterile, siliconized tubes and pipette tips.

II. Preparation of Stock Solutions

The positive control reagent contains 10⁴ copies/10 μ L of the HIV-1 control plasmid DNA and can be used directly. Further dilutions of the positive control can be prepared by diluting into either TE buffer (1 mM EDTA, 10 mM Tris-HCl, pH 8.0) or the negative control DNA.

If these control reagents are used with the GeneAmp PCR Reagent Kits with either native Taq DNA polymerase (Part No. N801-0043) or AmpliTaq[®] DNA polymerase (Part No. N801-0055), a sterile 10 mM MgCl₂ solution is needed, to supplement the MgCl₂ already present in the [10x] Reaction buffer. Dissolve 203 mg MgCl₂·6H₂O in 100 mL distilled water, filter through a 0.2 μ m filter and autoclave at 121°C at 15 PSI for 20 minutes.

III. Reaction Mix using the GeneAmp PCR Core Reagents (Part No. N808-0009)

Component	Addition Order	Volume	Final Concentration
Master Mix {	1	575 μ L	
	2	10 μ L	1X
	3	2 μ L	200 μ M
		2 μ L	200 μ M
		2 μ L	200 μ M
		2 μ L	200 μ M
	4	2 μ L	0.5 μ M
5	2 μ L	0.5 μ M	
6	0.5 μ L	2.5 Units/100 μ L	
7	10 μ L	2.5 mM	
Positive Control DNA		10 μ L	10 ⁴ copies/100 μ L
-or-			
Negative Control DNA	8	10 μ L	0.1 μ g human placental DNA/100 μ L
-or-			
Experimental sample		10 μ L*	
Total Mix		100 μL	

*Any combination of sterile distilled water and experimental sample volumes can be used as long as the total combined volume equals 67.5 μ L.

PERKIN ELMER CETUS

IV. Reaction Mix using the GeneAmp[®] PCR Reagent Kits (Part Nos. N801-0043 & N801-0055)

Component	Addition Order	Volume	Final Concentration
Master Mix	1	57.5 μ L	
	2	10 μ L	1X
	3	2 μ L	200 μ M
		2 μ L	200 μ M
		2 μ L	200 μ M
		2 μ L	200 μ M
	4	2 μ L	0.5 μ M
	5	2 μ L	0.5 μ M
	6	0.5 μ L	2.5 Units/100 μ L
	7	10 μ L	2.5 mM (including MgCl ₂ in [10x] reaction buffer)
Positive Control		10 μ L	10 ⁴ copies/100 μ L
-or-		-or-	
Negative Control	8	10 μ L	0.1 μ g human placental DNA/100 μ L
-or-		-or-	
Experimental sample		10 μ L	
Total Mix		100 μ L	

*Any combination of sterile distilled water and experimental sample volumes can be used as long as the total combined volume equals 67.5 μ L.

V. Notes on Reaction Mix for Controls

Vortex briefly (2 seconds; avoid generating bubbles) the *AmpliTag* or native *Taq* DNA Polymerase solutions or other recently thawed reagents, then spin down in a microcentrifuge before pipetting. Pipette enzyme carefully and slowly. The viscosity of the 50% glycerol in the buffer can lead to errors. If possible use a positive displacement pipette. Using master mixes (see above) will increase accuracy, reduce reagent loss on tips, and reduce tube-to-tube variability.

To reduce evaporation or refluxing, overlay the mix with with 50 to 100 μ L of mineral oil (Part No. 0186-2302). The oil should not interfere when withdrawing samples. If the entire volume is to be recovered, 100 μ L of high purity chloroform should be added after amplification. The aqueous phase containing the DNA will then float on the chloroform-oil mixture, allowing easy collection.

VI. Temperature Cycling for the Control Reagents

Optimum performance of the GeneAmp[™] PCR Process is achieved using the Perkin-Elmer Cetus DNA Thermal Cycler. For amplification using the GeneAmp[™] HIV-1 Reagents, program the Perkin-Elmer Cetus DNA Thermal Cycler for 4 linked files.

File # 11 STEP-CYCLE FILE PARAMETERS

Temp = 95C
 Seg #1 1min 0sec
 Temp = 0C
 Seg #2 0min 0sec
 Auto Segment Extension: OFF
 Cycle Count = 1
 Linked to file #12

File # 12 STEP-CYCLE FILE PARAMETERS

Temp = 95C
 Seg #1 1min 0sec
 Temp = 60C
 Seg #2 1min 0sec
 Temp = 0C
 Seg #3 0min 0sec
 Auto Segment Extension: OFF
 Cycle Count = 30
 Linked to file #13

File # 13 STEP-CYCLE FILE PARAMETERS

Temp = 60C
 Seg #1 10min 0sec
 Temp = 0C
 Seg #2 0min 0sec
 Auto Segment Extension: OFF
 Cycle Count = 1
 Linked to file #14

File # 14 SOAK FILE PARAMETERS

Temp = 4C

Refer to the Thermal Cycler User Manual for additional programming instructions.

The samples can be stored at 2° to 8°C until subsequent analysis. The selection of 60°C for the anneal-extend temperature is optimal for amplification of the positive control template. Due to genetic variation in the HIV-1 genome^{***} a lower anneal-extend temperature in the range of 52° to 60°C may be necessary to achieve amplification of viral sequences in experimental samples^{***}.

PERKIN ELMER CETUS

VII. Performance Claims

Using the amplification conditions described above and the Perkin-Elmer Cetus DNA Thermal Cycler and starting with 50 copies of the positive control template, sufficient amounts of the 115 bp PCR product will be generated to allow specific detection by the oligomer-hybridization procedure. Detailed procedures for labeling the SK19 probe with γ - 32 P]-ATP can be found in references 14 and 15.

VIII. Sequences

The sequences of the HIV-1 primers and probes and their location in the HIV-1 genome of isolate ARV-2¹¹. GenBank accession number KO2007, are shown below


Primer or Probe	Sequence (5'-3')	Location
SK38 Primer	ATAATCCACCTATCCCAGTAGGAGAAAT	<i>gag</i> 1551-1578
SK39 Primer	TTTGGTCCITGTCTTATGTCCAGAATGC	<i>gag</i> 1638-1665
SK19 Probe	ATCCTGGATTAAATAAAATAGTAAGAATGTATAGCCCTAC	<i>gag</i> 1595-1635

The North American consensus sequence of the 115 bp target¹² segment of the HIV-1 *gag* gene and the sequence of the 115 bp PCR product obtained using the SK38 and SK39 primers with the positive control plasmid DNA are shown below. Information on the extent of sequence conservation in this region among different HIV-1 isolates can be found in References 9 and 13.

	Sequence
North American Consensus	(5') *TAATCCACC TA*CCGAGTA GGAGAAAT*T ATAAAAG*TG GATAAT*CTG GG*TTAAATA AAATAGTAAG *ATGTATAG* CC* **CAGCA TTCT*GACAT AAGACAAGGA CC*AA (3') * nonconsensus nucleotide
115 bp PCR Product	< -----Region for SK38----- > (5') ATAATCCACC TATCCCAGTA GGAGAAATCT ATAAAAGATG GATAATCCTG GGATTAATA AAATAGTAAG AATGTATAGC CCTGTCAGCA TTCTGGACAT AAGACAAGGA CCAAA (3') < -----Region for SK39----- >

IX. HIV-1 Positive Control Plasmid

The following schematic diagram (not drawn to scale) shows the general rearrangement of the HIV-1 genome in the control plasmid. The *Bam* HI - *Nru* I region of plasmid pBR322 has been substituted with the disrupted (*pol* gene) and rearranged HIV26 sequence.



REFERENCES

- Ou, C-Y., Kwok, S., Mitchell, S.W., Mack, D.H., Sninsky, J.J., Krebs, J.W., Feorino, P., Warfield, D., Schochetman, G. 1988. DNA amplification for direct detection of HIV-1 in DNA of peripheral blood mononuclear cells. *Science* 239: 295-297.
- Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. 1990. *PCR Protocols: A Guide to Methods and Applications*. Academic Press, Inc., San Diego, CA.
- Mullis, K.B. and Faloona, F.A. 1987. Specific synthesis of DNA *in vitro* via a polymerase-catalyzed chain reaction. In: Wu, R. ed. *Methods in Enzymology*, vol. 155, Academic Press, San Diego, CA, pp. 335-350.
- Saiki, R.K., Scharf, S., Faloona, F., Mullis, K.B., Horn, G.T., Erlich, H.A., and Arnheim, N. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* 230: 1350-1354.
- Scharf, S.J., Horn, G.T., Erlich, H.A. 1986. Direct cloning and sequence analysis of enzymatically amplified genomic sequences. *Science* 233: 1076-1078.
- Saiki, R.K., Gelfand, D.H., Stoffel, S., Scharf, S.J., Higuchi, R., Horn, G.T., Mullis, K.B., and Erlich, H.A. 1988. Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase. *Science* 239: 487-491.
- Kawasaki, E.S., Clark, S.S., Coyne, M.Y., Smith, S.D., Champin, R., Whittle, O.N., and McCormick, F.P. 1988. Diagnosis of chronic myeloid and acute lymphocytic leukemias by detection of leukemia-specific mRNA sequences amplified *in vitro*. *Proc. Natl. Acad. Sci. USA*. 85: 5698-5702.

PERKIN ELMER CETUS

- 8 Higuchi, R. and Kwok, S. 1989. Avoiding false positives with PCR. *Nature* 339: 237-238.
- 9 Hahn, B.H., Gonda, M.A., Shaw, G.M., Popovic, M., Hoxie, J.A., Gallo, R.C. and Wong-Staal, F. 1985. Genomic diversity of the acquired immune deficiency syndrome virus HTLV-III: Different viruses exhibit greatest divergence in their envelope genes. *Proc. Natl. Acad. Sci. USA.* 82: 4813-4817.
- 10 Saag, M.S., Hahn, B.H., Gibbons, J., Li, Y., Parks, E.S., Parks, W.P. and Shaw, G.M. 1988. Extensive variation of human immunodeficiency virus type-1 *in vivo*. *Nature* 334: 440-444.
- 11 Goodenow, M., Huet, T., Saurin, W., Kwok, S., Sninsky, J. and Wain-Hobson, S. 1989. HIV-1 isolates are rapidly evolving quasispecies: Evidence for viral mixtures and preferred nucleotide substitutions. *Journal of Acquired Immune Deficiency Syndromes* 2: 344-352.
- 12 Meyernans, A., Cheyner, R., Albert, J., Seth, M., Kwok, S., Sninsky, J., Morfeldt-Manson, L., Asjo, B. and Wain-Hobson, S. 1989. Temporal fluctuations in HIV quasispecies *in vivo* are not reflected by sequential HIV isolations. *Cell* 58: 901-910.
13. Los Alamos Human Retroviruses and AIDS Directory. 1989. Los Alamos National Laboratory, Los Alamos, NM, 87545.
14. Kellogg, D. and Kwok, S. 1990. Detection of human immunodeficiency virus. In: Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. PCR Protocols: A Guide to Methods and Applications. Academic Press, Inc., San Diego, CA, pp. 337-347.
15. Molecular Cloning: A Laboratory Manual, 2nd Edition, Volume 2, Sambrook, J., Fritsch, E.F. and Maniatis, T., eds. 1989. Cold Spring Harbor Laboratory Press, N.Y., pp. 11.31-11.33.
16. Sanchez-Pescador, R., Power, M.D., Barr, P.J., Stermer, K.S., Stempien, M.M., Brown-Shimer, S.L., Gee, W.W., Renard, A., Randolph, A., Levy, J.A., Dina, D., and Luciw, P.A. 1985. Nucleotide sequence and expression of an AIDS-associated retrovirus (ARV2). *Science* 227: 484-492.
17. Hart, C., Chang, S-Y., Kwok, S., Sninsky, J., Ou, C.Y., and Schochetman, G. 1990. A replication-deficient HIV-1 DNA used for quantitation of the polymerase chain reaction (PCR). *Nucleic Acids Research* 18: 4029-4030.

PERKIN ELMER CETUS

ORDERING INFORMATION
PRODUCTS ARE NOT FOR DIAGNOSTIC OR MEDICAL USE
FOR RESEARCH USE ONLY

<p><i>In the United States:</i></p> <p>For PCR Technical support call toll-free 1-800-762-4001. FAX: 203-834-4391</p> <p>For easy ordering from PE XPRESS call toll-free 1-800-762-4002.</p> <p>To inquire about status of orders or to request a catalog from PE XPRESS call toll-free 1-800-762-4003.</p> <p>For product literature on the full line of Perkin-Elmer Cetus Instruments and Reagents call toll-free 1-800-762-4000.</p>	<p><i>Outside of the United States:</i></p> <p>For ordering information or the number of your local Perkin-Elmer sales representa- tive call:</p> <table style="width: 100%;"> <tr><td>Austria</td><td>(222) 623 101</td></tr> <tr><td>Belgium</td><td>02/21 20 611</td></tr> <tr><td>Canada</td><td>514-737-7575</td></tr> <tr><td>Denmark</td><td>042/811555</td></tr> <tr><td>Far East-Australia</td><td>3-560-4566</td></tr> <tr><td>Finland</td><td>90/880 144</td></tr> <tr><td>France</td><td>(1) 30 85 63 63</td></tr> <tr><td>Italy</td><td>39-83841</td></tr> <tr><td>Latin America-Mexico</td><td>5-651-7077</td></tr> <tr><td>Netherlands</td><td>0 1820-28122</td></tr> <tr><td>Norway</td><td>02/686565</td></tr> <tr><td>Portugal</td><td>1762172</td></tr> <tr><td>Spain</td><td>447-1021</td></tr> <tr><td>Sweden</td><td>08/7330010</td></tr> <tr><td>Switzerland</td><td>061/6925737</td></tr> <tr><td>United Kingdom</td><td>02407/4411</td></tr> <tr><td>West Germany</td><td>06103-708-118</td></tr> </table>	Austria	(222) 623 101	Belgium	02/21 20 611	Canada	514-737-7575	Denmark	042/811555	Far East-Australia	3-560-4566	Finland	90/880 144	France	(1) 30 85 63 63	Italy	39-83841	Latin America-Mexico	5-651-7077	Netherlands	0 1820-28122	Norway	02/686565	Portugal	1762172	Spain	447-1021	Sweden	08/7330010	Switzerland	061/6925737	United Kingdom	02407/4411	West Germany	06103-708-118
Austria	(222) 623 101																																		
Belgium	02/21 20 611																																		
Canada	514-737-7575																																		
Denmark	042/811555																																		
Far East-Australia	3-560-4566																																		
Finland	90/880 144																																		
France	(1) 30 85 63 63																																		
Italy	39-83841																																		
Latin America-Mexico	5-651-7077																																		
Netherlands	0 1820-28122																																		
Norway	02/686565																																		
Portugal	1762172																																		
Spain	447-1021																																		
Sweden	08/7330010																																		
Switzerland	061/6925737																																		
United Kingdom	02407/4411																																		
West Germany	06103-708-118																																		
<p><i>In Europe:</i></p> <p>For PCR Technical support call Perkin-Elmer's European Sales Support Center in Vaterstetten, West Germany at 49/8106/381-112. FAX: 49/8106/6697</p>																																			

<p>N801-0150 DNA Thermal Cycler, 120 volts, 60 Hz N801-0177 DNA Thermal Cycler, 240 volts, 50 Hz N801-0189 DNA Thermal Cycler, 100 volts, 50/60 Hz</p>	<p>} A microprocessor controlled, thermal cycling instrument which automates the rapid and precise temperature changes needed in the GeneAmp™ PCR process.</p>
--	--

Perkin-Elmer® is a registered trademark of Perkin-Elmer Corporation.
Cetus® is a registered trademark of Cetus Corporation.
GeneAmp® and AmpliTaq® are registered trademarks of Cetus Corporation, licensed to Perkin-Elmer Cetus Instruments.
GeneAmplicon™ is a trademark of Cetus Corporation, licensed exclusively to Perkin-Elmer Cetus Instruments.


Copyright 1990 by Perkin-Elmer Cetus Instruments.

BIO-24
55905 9/90

11/11/91

Four copies of the book sent via messenger to David McBride on 11/10/91. David will have these copies inserted into final Exhibit Package for Court, Berger, Jones, Day.

Ross



ATTORNEY'S EYES ONLY
RESTRICTED

RMS 71157

PERKIN ELMER CETUS
761 Main Ave., Norwalk, CT 06859

GeneAmp® PCR Carry-over Prevention Kit

NOT FOR MEDICAL OR DIAGNOSTIC USE
FOR RESEARCH USE ONLY
See notice to purchaser

Part No. N808-0068

The capability for single molecule amplification provided by the GeneAmp™ Polymerase Chain Reaction (PCR) process^{1,2} has required the adoption of special laboratory practices to avoid false positive amplifications.³ A primary source of exogenous templates leading to false positive amplifications has been identified as products from previous PCR amplifications.⁴ The Perkin-Elmer Cetus GeneAmp® Carry-over Prevention Kit has been developed to provide a means for ensuring that products of previous PCR amplifications cannot be reamplified to produce false positive results. The GeneAmp Carry-over Prevention Kit uses enzymatic and chemical reactions analogous to the restriction-modification and excision-repair systems of cells to degrade specifically PCR products from previous PCR amplifications, but not degrade native nucleic acid templates. The method used to make PCR products susceptible to degradation involves substituting dUTP for dTTP in the PCR reaction mix and treating subsequent PCR amplifications with the enzyme Uracil N-glycosylase (UNG, EC 3.2.2.-) prior to amplification.

The UNG provided in this kit is a pure, nuclease free, 26 kDa enzyme encoded by the *Escherichia coli* Uracil N-glycosylase gene⁵ which has been inserted into an *E. coli* host to direct the expression of the native form of the enzyme.⁶ When dUTP replaces dTTP as a dNTP substrate in PCR and the method described below is used, UNG treatment can prevent the subsequent reamplification of dU-containing PCR products.⁷ This insert contains important details about the use of dUTP and UNG and the properties of dU-containing oligonucleotides. *The entire product insert should be read prior to using the reagents provided.*

INTRODUCTION

Bacterial cells are particularly efficient at (i) degrading incoming phage DNA if the viral genomes are not methylated^{8,9} and (ii) removing by enzymatic excision modified bases that result from inappropriate insertion by the replication apparatus or from damage by mutagens.¹⁰ Warner and colleagues, for example, demonstrated that T4¹¹ and T5¹² phage containing about 30% of the thymine replaced by uracil were restricted by wild type *E. coli* but were capable of nearly equivalent infection as phage lacking uracil in *E. coli* cells deficient in UNG. In addition, Kunkel¹³ exploited the *in vivo* degradation of uracil-containing DNA to enhance the recovery of recombinant M13 phage from *in vitro* mutagenesis reactions. The GeneAmp PCR Carry-over Prevention Kit adapts the UNG and dU-containing DNA cellular restriction-excision system to provide a means of degrading PCR products from previous reactions.




Figure 1. Cleavage of dU-containing PCR products by UNG in concert with alkaline pH and heat.

PERKIN ELMER CETUS

The distinguishing characteristic of PCR product relative to native DNA template is that PCR product is newly synthesized and thus modified nucleotides other than the typical four nucleotides (dA, dG, dT and dC) can be incorporated. Incorporation of a novel or modified nucleotide in the PCR product thus serves as the basis for distinguishing the newly synthesized PCR product from the starting native DNA template. Further, since enzymes of the excision-repair system can either remove the altered base and/or cleave the polynucleotide backbone, specific degradation of the PCR product instead of the starting native DNA template can be achieved. Abasic polynucleotides are susceptible to hydrolysis in alkaline solutions especially at elevated temperatures.¹ The alkaline pH of the PCR buffer and the elevated temperature of the initial denaturation step (> 94°C) are more than sufficient to result in abasic polynucleotide strand scission. The reagents in this kit can be used to modify PCR amplifications such that the PCR products from previous amplifications are susceptible to enzymatic excision and chemical cleavage without degrading native DNA templates.

To render selectively PCR products susceptible to cleavage, the modification of PCR amplifications begins with the substitution of dUTP for dTTP in the PCR reaction mix. UNG is added to the PCR reaction mix. A short room temperature incubation step is performed prior to the start of the PCR. This allows time for UNG to excise the uracil in carried-over dU-containing PCR products should they be present in the reaction. The excised dU-containing PCR products are refractory to further PCR amplification because of the frequent stalling of DNA polymerases at the abasic sites^{2,3} and/or cleavage of the abasic polynucleotide strands due to the inherent alkaline and thermal lability of the resulting abasic linkage. Since UNG is active on single- and double-stranded DNA⁴, the procedure should work on dU-containing PCR products from standard or asymmetric PCR amplifications. Ribouracil residues in RNA and dUTP are not substrates for UNG.⁵ Figure 1 shows the reactions involved in the UNG catalyzed degradation of dU-containing DNA.

LIST OF COMPONENTS

Store UNG and dUTP solutions at -20°C in a constant temperature freezer. If stored under the proper conditions, the product will maintain performance until the control date printed on the label.

Reagent	Volume	Concentration	Comments
Uracil N-glycosylase	100 µL	1 Unit/µL	See Uracil N-glycosylase Characteristics
dUTP	320 µL	20 mM	Deoxynucleotide dissolved in glass distilled water; titrated with NaOH to pH 7.0

NOTICE TO PURCHASER: DISCLAIMER OF LICENSE

This product is optimized for use in the *GeneAmp*[™] Polymerase Chain Reaction (PCR) process covered by patents issued to Cetus Corporation, but no license under these patents to use the *GeneAmp* PCR process is conveyed by the purchase of this product. A license to use the *GeneAmp* PCR process for certain research and development activities accompanies the purchase and use of Perkin-Elmer Cetus' *GeneAmp* PCR Core Reagents, *GeneAmp* PCR Reagent Kits, and certain other products, or is available from Cetus Corporation. For information contact Cetus Director of Business Development, PCR Division, at 1400 Fifty-Third Street, Emeryville, CA 94608. Tel. (415) 420-3300.

URACIL N-GLYCOSYLASE CHARACTERISTICS

Concentration:	1 Unit/µL
Unit Definition:	One unit of enzyme is defined as the amount that releases 1 nmole of uracil from a dU-containing DNA template into acid soluble material per 60 minutes at 37°C under the analysis conditions below.
UNG Activity Analysis Conditions:	10 mM Tris-HCl, pH 8.3 (at room temperature); 50 mM KCl; and 100 pmoles dUMP substrate (dU-containing Lambda PCR product, a mix of unlabeled and [³ H-Uracil]-labeled), in a final volume of 50 µL and incubated at 37°C for 10 min. The dU-containing 500 bp Lambda PCR product is prepared by following the procedure for amplifying the Control Template provided with the <i>GeneAmp</i> [®] PCR Reagent Kit (Part No. N801-0055) except that a mixture of dUTP and [³ H-Uracil]-dUTP are substituted for dTTP.
Storage Buffer:	5% Glycerol (v/v) 150 mM NaCl 30 mM Tris-HCl, pH 7.5 (at room temperature) . . 1.0 mM EDTA (Ethylenediamine-tetraacetic acid) 1 mM DTT (Dithiothreitol) 0.05% Tween 20 [®]
Storage Temperature:	Store UNG at -20°C, in a constant temperature freezer. If stored under proper conditions, the enzyme will remain active until the control date printed on the label.
Associated Activities:	No single-stranded or double-stranded endonuclease activities were detected when 5 Units of UNG were incubated at 37°C for 1 hour with 1.2 µg of circular M13 DNA or 600 ng supercoiled pBR322 (dcm ⁺ , dam ⁺), respectively, in a final volume of 50 µL. No double-stranded exonuclease activity was detected when 5 Units of UNG were incubated at 37°C for 30 minutes with 1 pmole of [³ H]-dT-labeled 500 bp Lambda PCR product in a final volume of 50 µL. No 5'-single-stranded exonuclease activity was detected when 5 Units of UNG were incubated at 37°C for 30 minutes with 0.3 pmoles of a 5'-[³² P]-labeled oligonucleotide (41-mer) in a final volume of 50 µL. No 3'-single-stranded exonuclease activity was detected when 5 Units of UNG were incubated at 37°C for 30 minutes with an oligonucleotide (40-mer) labeled at the 3' end with α-[³² P]-Cordycepin in a final volume of 50 µL.

PERKIN ELMER CETUS

PROTOCOL FOR DNA AMPLIFICATION

I. Materials Necessary but Not Provided

The following materials are needed in order to amplify target DNA segments with dUTP

- *GeneAmp*[®] PCR Core Reagents (Part No. N808-0009)
- or-
- *GeneAmp* PCR Reagent Kit with *AmpliTaq*[®] DNA Polymerase (Part No. N801-0055)
- or-
- *GeneAmp* PCR Reagent Kit with Native *Taq* DNA Polymerase (Part No. N801-0043)
- and-
- Primers and Experimental template.

II. General Advice

Because of the enormous amplification possible with the *GeneAmp*[™] PCR process, small levels of DNA carry over from samples with high DNA levels, positive control templates, or from previous PCR amplifications can result in product even in the absence of purposefully added template DNA. Although the protocol and reagents described here are capable of degrading or eliminating large numbers of carried over PCR products, we encourage users to continue using the specific devices and suggestions described in the package inserts provided with the products listed in Section I and in References 4 and 5 to minimize cross-contamination from non-dU-containing PCR products or other samples. Briefly, the use of dedicated or disposable vessels, solutions, and pipettes (preferably positive displacement pipettes) for DNA preparation, reaction mixing, and sample analysis is still recommended. In addition, all reactions should be set up in an area separate from PCR product analysis.

A master mix of reagents (water, buffer, dNTPs [dATP, dCTP, dGTP, and dUTP], primers, *AmpliTaq*[®] DNA Polymerase, and UNG) for all samples can be prepared first, then aliquoted to individual tubes. Magnesium chloride and template DNA are then added. Using such mixes will minimize pipetting losses, increase accuracy, and reduce the number of reagent transfers.

Amplifications are performed in 100 µL of reaction mix in capped 0.5 mL polypropylene microcentrifuge tubes. Perkin-Elmer Cetus 0.5 mL *GeneAmp*[™] Reaction Tubes (Part No. N801-0180) provide the best heat transfer when using the Perkin-Elmer Cetus DNA Thermal Cycler, because of their uniform fit in the wells. Since DNA may stick to plastic and since nucleases are often found on surfaces, it may be preferable to use sterile, siliconized tubes and pipette tips.

III. General Protocol using the *GeneAmp* PCR Carry-over Prevention Kit and *GeneAmp* PCR Core Reagents

Reaction Mix

Component	Addition Order	Volume	Final Concentration
Sterile Distilled Water	1	-	-
10X PCR Buffer II	2	10 µL	1X
dATP, 10 mM	3	2 µL	200 µM
dCTP, 10 mM		2 µL	200 µM
dGTP, 10 mM		2 µL	200 µM
dUTP, 20 mM		1.5 µL	200-1000 µM**
Primer 1	4	1.5 µL	0.2-1.0 µM
Primer 2	5	1.5 µL	0.2-1.0 µM
<i>AmpliTaq</i> [®] DNA Polymerase	6	0.5 µL	2.5 Units/100 µL
UNG	7	1 µL	1 Unit/100 µL***
25 mM MgCl ₂	8	4-16 µL	1.0-4.0 mM****
Experimental template	9	-	<1 µg/100 µL
Total Mix		100 µL	

- * Any combination of sterile distilled water and experimental template volumes can be used as long as the total volume of the reaction (including buffer, dNTPs, primers, enzymes and MgCl₂ solutions) equals 100 µL.
- ** Although many PCR systems are amplified efficiently when an equimolar concentration of dUTP (200 µM) is substituted for dTTP (200 µM) in the standard reaction mix, some PCR systems require increasing the dUTP concentration 2 to 5 fold (e. g., 400-1000 µM) for optimal amplification efficiencies. Increase MgCl₂ on an equimolar basis, when an increased dUTP concentration is required. In some systems, the level of amplification obtained even with 1000 µM dUTP may be less than the level obtained with 200 µM dTTP. Sufficient 20 mM dUTP solution is provided for 100 PCR amplifications (final volume of 100 µL) at 600 µM dUTP.
- *** Depending on the level of potential dU-containing PCR product cross-contamination, the required amount of UNG may vary. With a short incubation step (0-10 minutes at room temperature) 1 Unit/100 µL reaction is generally more than sufficient to preclude subsequent reamplification of high levels (1 x 10⁶ copies) of dU-containing PCR product. If further optimization of the concentration of UNG is desired, concentrations in the range of 1.0 to 2.0 Units/100 µL have been shown to be effective (see Section IV.3).
- **** The optimal magnesium chloride concentration may vary, depending on the total dNTP concentration and on the primer and templates used and must be determined empirically. In most cases a final concentration of MgCl₂ in the range of 1.0 to 5.0 mM in the reaction mix will work well. The addition of UNG will not affect the magnesium chloride optimum of your specific PCR amplification. Magnesium ion is not required for UNG activity.

UNG Reaction Conditions

After the addition of MgCl₂ and Experimental template DNA begin with a 1 to 10 minute incubation at room temperature (18° to 22°C) to allow UNG to excise uracil from any contaminating dU-containing PCR product from previous amplifications. The time of this incubation depends on the level of UNG used and the potential for primer-dimer formation in your the PCR system (see Section IV.3). Follow the room temperature incubation with at least a 10 minute incubation at 95°C using a soak file in the DNA Thermal Cycler process. Incubation at 95°C for at least 10 minutes is necessary to cleave the abasic dU-containing PCR product generated in the room temperature incubation, to inactivate UNG and to denature the native DNA in the experimental sample. If a

PERKIN ELMER CETUS

shorter incubation at 95°C is used. UNG may not be completely inactivated, leading to degradation of newly synthesized dU-containing PCR product during the PCR process. The 10 minute incubation at 95°C does not affect significantly the activity of *AmpliTag²* DNA Polymerase, which has a half-life of >45 minutes at 95°C.

Since UNG has activity below 55°C, the annealing temperature used for the PCR amplification should be at or above this temperature, in order to avoid degradation of newly synthesized dU-containing PCR products by residual UNG activity.

VERY IMPORTANT

After the final PCR cycle, use a soak file to hold PCR reaction mixes in which dUTP and UNG are present at 72°C until they are removed from the Perkin-Elmer Cetus DNA Thermal Cycler. Store reaction mixes immediately at -20°C or add an equal volume of chloroform to prevent any degradation of dU-containing PCR products by residual or reactivated UNG.


IV. Notes on Preventing PCR Product Carry-over

1. One should determine if the sequence being amplified contains dA and dT nucleotides. Only DNA sequences that have dA and dT will lead to dU-containing PCR products that can be degraded.
2. The primers that are used should contain dA nucleotides near the 3' termini so that any primer-dimer generated is efficiently degraded by UNG at least as well as any dU-containing PCR products. The further a dA nucleotide is from the 3' terminus the more likely that partially degraded primer-dimer molecules may serve as templates for a subsequent PCR amplification. Production of primer-dimer could compromise the amplification of the desired target region. If primers can not be selected with dA nucleotides near the termini, the use of primers with 3' terminal dU nucleotides should be considered. Terminal dU nucleotides are not substrates for UNG¹¹ and thus the primers will not be degraded. Biotin-dUMP derivatives are not substrates for UNG.
3. The incubation time of the UNG reaction may affect the subsequent PCR amplification efficiency. When 1 Unit/100 µL of UNG is used, a 10 minute incubation at room temperature will be sufficient to degrade 10⁸ copies of cross-contaminating dU-containing PCR products. *The concentration of UNG and the time of the incubation step necessary to prevent amplification of contaminating dU-containing PCR products depends on the PCR conditions necessary to amplify your particular DNA sequence and the level of contamination to be expected.* In most cases, the general recommendation of using UNG at 1 Unit/100 µL reaction and incubation at room temperature for 10 minutes will suffice. If the PCR system is prone to primer-dimer formation, a 10 minute incubation at room temperature may result in the preferential amplification of primer-dimer sequences instead of the desired target sequence. If primer-dimer amplification dominates the PCR, reduce the incubation time at room temperature with or without increasing the UNG concentration to 2 Units/100 µL reaction.

V. Temperature Cycling and Cycle Optimization

The performance of the *GeneAmp[®]* PCR Carry-over Prevention Kit is optimized and quality control tested for performance on the Perkin-Elmer Cetus DNA Thermal Cycler. Optimal performance of the *GeneAmp[®]* PCR process is influenced by choice of temperature, time at temperature, and length of time between temperatures for each step in the cycle.

A typical cycling profile for utilizing UNG to degrade contaminating dU-containing PCR products from previous amplifications prior to PCR amplification of freshly added native DNA Experimental Templates is shown below:

**VERY IMPORTANT**

After the final PCR cycle, use a soak file to hold PCR reaction mixes in which dUTP and UNG are present at 72°C until they are removed from the Perkin-Elmer Cetus DNA Thermal Cycler. Store reaction mixes immediately at -20°C or add an equal volume of chloroform to prevent any degradation of dU-containing PCR product by residual or reactivated UNG.

VI. Properties of dU-containing DNA

The substitution of dUTP for dTTP in PCR amplifications leads to the generation of PCR products containing dU instead of dT. Although this renders the dU-containing PCR products susceptible to degradation by UNG, most other properties of dU-containing PCR products are largely not affected.

- dU-containing PCR products will serve in an equivalent manner as dT-containing PCR products as hybridization targets.¹²⁻¹⁴
- dU-containing PCR products should also serve as templates for dideoxy-terminated sequencing reactions.
- dU-containing DNA can be cloned directly, if transformed into a *ung⁻* bacterial host.^{15, 16}
- The recognition of dU-containing DNA by restriction endonucleases has been studied.¹⁷⁻²⁰ Depending on the specific endonuclease there may be no effect of the substitution of dU for dT on enzyme activity (e.g., *EcoRI* and *BamHI*) or the dU-containing DNA is cleaved more slowly than dT-containing DNA (e.g., *HpaI*, *HindII*, and *HindIII*). For other endonucleases the effect of substituting dU for dT on enzyme activity will need to be examined empirically on an individual enzyme basis.

The use of PCR products containing exclusively uracil is not recommended for protein binding or site-specific recognition studies (e.g., operators²¹ and promoters²², restriction endonucleases^{23, 24}, or methylases). The absence of the C5 methyl group of thymine in the major groove, the face of the helix most frequently involved in protein-nucleic acid interactions^{25, 26}, and/or distortion of the glycosidic bonds of the dA:dU (and adjacent base pairs)²⁷ may impair recognition by proteins.

PERKIN ELMER CETUS

REFERENCES

1. Mullis, K.B., and Faloona, F.A. 1987. Specific synthesis of DNA *in vitro* via a polymerase-catalyzed chain reaction. In: Wu, R. ed. *Methods in Enzymology*, Vol 155, Academic Press, San Diego, CA, pp. 335-350.
2. Saiki, R.K., Scharf, S.J., Faloona, F., Mullis, K.B., Horn, G.T., Erlich, H.A., and Arnheim, N. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* 230: 1350-1354.
3. Saiki, R.K., Gelfand, D.H., Stoffel, S., Scharf, S.J., Higuchi, R., Horn, G.T., Mullis, K.B., and Erlich, H.A. 1988. Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase. *Science* 239: 487-491.
4. Higuchi, R. and Kwok, S. 1989. Avoiding false positives with PCR. *Nature* 339: 237-238.
5. Kwok, S. 1990. Chapter 17. Procedures to minimize PCR-product carry-over and Orrego, C. 1990. Chapter 54. Organizing a laboratory for PCR work. In: Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. *PCR Protocols. A Guide to Methods and Applications*. Academic Press, Inc., San Diego, CA, pp. 142-145 and 447-454, respectively.
6. Varshney, U., Hutcheon, T., and van de Sande, J.H. 1988. Sequence analysis, expression, and conservation of *Escherichia coli* uracil DNA glycosylase and its gene (*ung*). *J. Biol. Chem.* 263: 7776-7784.
7. Sninsky, J.J., Gates, C., McKinney, N., Birch, D., Akers, J., Lawyer, F.C., and Gelfand, D.H. 1990. Manuscript submitted.
8. Bertani, G., and Weigle, J.J. 1953. Host controlled variation in bacterial viruses. *J. Bacteriol.* 65: 113-121.
9. Arber, W., and Dussoix, D. 1962. Host specificity of DNA produced by *Escherichia coli*. I. Host controlled modification of bacteriophage λ . *J. Mol. Biol.* 5: 18-36.
10. Dussoix, D., and Arber, W. 1962. Host specificity of DNA produced by *Escherichia coli*. II. Control over acceptance of DNA from infecting phage λ . *J. Mol. Biol.* 5: 37-49.
11. Lindahl, T. 1976. New class of enzymes acting on damaged DNA. *Nature* 259: 64-66.
12. Warner, H.R. and Duncan, B.K. 1978. *In vivo* synthesis and properties of uracil-containing DNA. *Nature* 272: 32-34
13. Warner, H.R., Thompson, R.B., Mozer, T.J., and Duncan, B.K. 1979. The properties of a bacteriophage T5 mutant unable to induce deoxyuridine 5'-triphosphate nucleotidohydrolase: synthesis of uracil-containing T5 deoxyribonucleic acid. *J. Biol. Chem.* 254: 7534-7539.
14. Kunkel, T. 1985. Rapid and efficient site-specific mutagenesis without phenotypic selection. *Proc. Natl. Acad. Sci. USA* 82: 488-492.
15. Tamm C., Shapiro, H.S., Lipshitz, R., and Chargaff, E. 1953. Distribution density of nucleotides within a desoxyribonucleic acid chain. *J. Biol. Chem.* 203: 673-688.
16. Lindahl, T., and Andersson, A. 1972. Rate of chain breakage at apurinic sites in double-stranded deoxyribonucleic acid. *Biochemistry* 11: 3618-3623.
17. Sagher, D. and Strauss, B. 1983. Insertion of nucleotides opposite apurinic/apyrimidinic sites in deoxyribonucleic acid during *in vitro* synthesis: uniqueness of adenine nucleotides. *Biochemistry* 22: 4518-4526.
18. Takeshita, M., Chang, C-N., Johnson, F., Will, S., and Grollman, A.P. 1987. Oligodeoxynucleotides containing synthetic abasic sites. Model substrates for DNA polymerases and apurinic/apyrimidinic endonucleases. *J. Biol. Chem.* 262: 10171-10179.
19. Schaaper, R.M., Kunkel, T.A., and Loeb, L.A. 1983. Infidelity of DNA synthesis associated with bypass of apurinic sites. *Proc. Natl. Acad. Sci. USA* 80: 487-491.
20. Randall, S.K., Eritja, R., Kaplan, B.E., Patruska, J., and Goodman, M.F. 1987. Nucleotide insertion kinetics opposite abasic lesions in DNA. *J. Biol. Chem.* 262: 6864-6870.
21. Lindahl, T., Ljungquist, S., Siegart, W., Nyberg, B., and Sperens, B. 1977. DNA glycosidases. Properties of uracil-DNA glycosidase from *Escherichia coli*. *J. Biol. Chem.* 252: 3286-3294.
22. Delort, A-M., Duplax, A-M., Molko, D., and Teoule, R. 1985. Excision of uracil residues in DNA: mechanism of action of *Escherichia coli* and *Micrococcus luteus* uracil-DNA glycosylases. *Nuc. Acids Res.* 13: 319-335.
23. Brennan, C.A., and Gumport, R.I. 1985. T4 RNA ligase catalyzed synthesis of base analogue-containing oligodeoxyribonucleotides and a characterization of their thermal stabilities. *Nuc. Acids Res.* 24: 8665-8684.
24. Berkner, K.L., and Folk, W.R. 1979. The effects of substituted pyrimidines in DNAs on cleavage by sequence-specific endonucleases. *J. Biol. Chem.* 254: 2551-2560.
25. Bodnar, J.W., Zempsky, W., Warder, D., Bergson, C., and Ward, D.C. 1983. Effect of nucleotide analogs on the cleavage of DNA by the restriction enzymes *AclI*, *DdeI*, *HinfI*, *RsaI*, and *TaqI*. *J. Biol. Chem.* 258: 15206-15213.
26. Marchionni, M.A., and Roufa, D.J. 1978. Digestion of 5-Bromodeoxyuridine-substituted λ -DNA by restriction endonucleases. *J. Biol. Chem.* 253: 9075-9081.
27. Goeddel, D.V., Yansura, D.G., and Caruthers, M.H. 1977. Studies on gene control regions. VI. The 5-methyl of thymine, a lac repressor recognition site. *Nuc. Acids Res.* 9: 3039-3054.
28. Dubendorff, J.W., deHaseth, P.L., Rosendahl, M.S., and Caruthers, M.H. 1987. DNA functional groups required for formation of open complexes between *Escherichia coli* RNA polymerase and the λ P_q promoter. Identification via base analog substitutions. *J. Biol. Chem.* 262: 892-898.
29. Ivarie, R. 1987. Thymine methyls and DNA - protein interactions. *Nuc. Acids Res.* 15: 9975- 9983.
30. Seeman, N.C., Rosenberg, J.M., and Rich, A. 1976. Sequence-specific recognition of double helical nucleic acids by proteins. *Proc. Natl. Acad. Sci. USA* 73: 804-808.
31. Delort, A-M., Neumann, J.M., Molko, D., Herve, M., Teoule, R., and Dinh, S.T. 1985. Influence of uracil defect on DNA structure: ¹H NMR investigation at 500 MHz. *Nuc. Acids Res.* 9: 3343- 3355.
32. Duncan, B.K. 1985. Isolation of insertion, deletion, and nonsense mutations of the uracil-DNA glycosylase (*ung*) gene of *Escherichia coli* K-12. *J. Bacteriol.* 164: 689-695.

PERKIN ELMER CETUS

ORDERING INFORMATION

PRODUCTS ARE NOT FOR MEDICAL OR DIAGNOSTIC USE
FOR RESEARCH USE ONLY

<p>In the United States:</p> <p>For PCR Technical support call toll-free 1-800-762-4001. FAX: 203-834-4391</p> <p>For easy ordering from PE XPRESS call toll-free 1-800-762-4002.</p> <p>To inquire about status of orders or to request a catalog from PE XPRESS call toll-free 1-800-762-4003.</p> <p>For product literature on the full line of Perkin-Elmer Cetus Instruments and Reagents call toll-free 1-800-762-4000.</p>	<p>Outside of the United States:</p> <p>For ordering information or the number of your local Perkin-Elmer sales representa- tive call:</p> <table> <tr><td>Austria</td><td>(222) 623 101</td></tr> <tr><td>Belgium</td><td>02/21 20 611</td></tr> <tr><td>Canada</td><td>514-737-7575</td></tr> <tr><td>Denmark</td><td>042/811555</td></tr> <tr><td>Far East-Australia</td><td>3-560-4566</td></tr> <tr><td>Finland</td><td>90/880 144</td></tr> <tr><td>France</td><td>(1) 30 85 63 63</td></tr> <tr><td>Italy</td><td>39-83841</td></tr> <tr><td>Latin America-Mexico</td><td>5-651-7077</td></tr> <tr><td>Netherlands</td><td>0 1820-28122</td></tr> <tr><td>Norway</td><td>02/686565</td></tr> <tr><td>Portugal</td><td>1762172</td></tr> <tr><td>Spain</td><td>447-1021</td></tr> <tr><td>Sweden</td><td>08/7330010</td></tr> <tr><td>Switzerland</td><td>061/6925737</td></tr> <tr><td>United Kingdom</td><td>02407/4411</td></tr> <tr><td>Germany</td><td>06103-708-118</td></tr> </table>	Austria	(222) 623 101	Belgium	02/21 20 611	Canada	514-737-7575	Denmark	042/811555	Far East-Australia	3-560-4566	Finland	90/880 144	France	(1) 30 85 63 63	Italy	39-83841	Latin America-Mexico	5-651-7077	Netherlands	0 1820-28122	Norway	02/686565	Portugal	1762172	Spain	447-1021	Sweden	08/7330010	Switzerland	061/6925737	United Kingdom	02407/4411	Germany	06103-708-118
Austria	(222) 623 101																																		
Belgium	02/21 20 611																																		
Canada	514-737-7575																																		
Denmark	042/811555																																		
Far East-Australia	3-560-4566																																		
Finland	90/880 144																																		
France	(1) 30 85 63 63																																		
Italy	39-83841																																		
Latin America-Mexico	5-651-7077																																		
Netherlands	0 1820-28122																																		
Norway	02/686565																																		
Portugal	1762172																																		
Spain	447-1021																																		
Sweden	08/7330010																																		
Switzerland	061/6925737																																		
United Kingdom	02407/4411																																		
Germany	06103-708-118																																		
<p>In Europe:</p> <p>For PCR Technical support call Perkin-Elmer's European Sales Support Center in Vaterstetten, Germany at 49/8106/381-112. FAX: 49/8106/6697</p>																																			

RELATED PRODUCTS FROM Perkin-Elmer Cetus:

Part Number	Description
N801-0150	DNA Thermal Cycler, 120 V, 60 Hz A microprocessor controlled thermal cycling instrument which automates the rapid and precise temperature changes needed in the GeneAmp PCR process.
N801-0177	DNA Thermal Cycler, 240 V, 50 Hz
N801-0189	DNA Thermal Cycler, 100 V, 50/60 Hz

NEW INSTRUMENTS FROM Perkin-Elmer Cetus:

N801-0100	DNA Thermal Cycler 480, 120 V, 60Hz
N801-0101	DNA Thermal Cycler 480, 240 V, 50 Hz
N801-0102	DNA Thermal Cycler 480, 100 V, 50/60 HZ

Perkin-Elmer® is a registered trademark of The Perkin-Elmer Corporation.
Cetus® is a registered trademark of the Cetus Corporation.
GeneAmp®, Ampflaq® and GeneAmp® are registered trademarks/trademark of Cetus Corporation, licensed to Perkin-Elmer Cetus Instruments.
Tween 20® is a registered trademark of ICI.
Copyright 1990 by Perkin-Elmer Cetus Instruments. All rights reserved.

BIO-42
55605-11/90

ATTORNEY'S EYES ONLY
RESTRICTED

RMS 71163