

GUIDANT

August 17, 2001

PMA Document Mail Center (HFZ-401)
Center for Devices and Radiological Health
Food and Drug Administration
9200 Corporate Boulevard
Rockville, MD 20850

Re: Annual Report for PMA P960040 - VENTAK® AV™ and VENTAK PRIZM
DR/VR Automatic Implantable Cardioverter Defibrillator

Enclosed is the annual report for the VENTAK AV family of Automatic Implantable Cardioverter Defibrillators (AICD™) as required in the Conditions of Approval Letter for PMA P960040 (dated July 18, 1997) and subsequent supplements to the PMA. The information in this report covers the period between June 1, 2000 and May 31, 2001.

GUIDANT considers the information contained in this report to be confidential information under 520(c) of the Act, and not releasable to the public domain for any reason.

If you have questions regarding this report, please contact me at 651-582-5224.

Sincerely,

GUIDANT CPI

[
Vanessa D. Ware
Regulatory Affairs Associate
GUIDANT - Cardiac Pacemakers (CPI)

RECEIVED

AUG 20 11 27 AM '01

FDA/CDRH/CDE/DHC

POST-APPROVAL ANNUAL REPORT
For
VENTAK AV and VENTAK PRIZM DR/VR
AUTOMATIC IMPLANTABLE CARDIOVERTER DEFIBRILLATOR SYSTEM

PMA NO P960040

GUIDANT - Cardiac Pacemakers (CPI)
4100 Hamline Avenue North
St. Paul, MN 55112

Submission Date: August 17, 2001

Prepared and Submitted by:

[Signature]

Vanessa D. Ware
Regulatory Affairs Associate

[Signature]

Nancy L. Formosa
Regulatory Affairs Associate

Approved by:

[Signature] For Gayle Johnson

Gayle Johnson
Manager, Regulatory Affairs

TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
I. REPORT OVERVIEW	2
A. Scope of Report	2
B. Recall and Safety Alert Update	3
II. DEVICE CHANGES	4
A. PMA Supplemental Submissions	4
B. Changes Not Included in a Submission	6
III. BIBLIOGRAPHY AND SUMMARY OF STUDIES	6
A. Unpublished Reports Regarding the Device	7
B. Reports in Scientific Literature Concerning the Device	8
IV. SYSTEM EXPERIENCE	22
V. LIFE TESTING RESULTS - VENTAK VR	23
APPENDICES:	
APPENDIX A	24
APPENDIX B	25
APPENDIX C	26

I. REPORT OVERVIEW

A. Scope of Report

This Post-Approval Annual Report presents performance information for the VENTAK® AV™ and VENTAK PRIZM DR/VR Automatic Implantable Cardioverter Defibrillator System as required in the Conditions of Approval Letter for PMA P960040, (Appendix A) and subsequent supplements.

Per the requirements of the Condition of Approval Letter, this annual report includes information on device performance for the period between June 1, 2000 and May 31, 2001. The cumulative survival analysis is an ongoing accumulation of data on device life since the first documented implant. The cumulative survival analysis in this report is based on the data collected since the first documented domestic implant through the period ending May 31, 2001.

The following models are included in this report:

Guidant Brand Name	Model #	Application Software
VENTAK AV AICD Pulse Generator	1810, 1815	2833
VENTAK AV II DDD AICD Pulse Generator	1820, 1825	2833
VENTAK AV II DR AICD Pulse Generator	1821, 1826	2843
VENTAK AV III DR AICD Pulse Generator	1831, 1836	2843
VENTAK VR AICD Pulse Generator	1774, 1775	2841
VENTAK PRIZM DR/VR AICD Pulse Generator	1850, 1851, 1855, 1856	2844
VENTAK PRIZM DR/HE AICD Pulse Generator	1853, 1858	2844
VENTAK PRIZM VR/HE AICD Pulse Generator	1852, 1857	2844
VENTAK PRIZM 2 DR AICD Pulse Generator	1861, 1862	2844

We have reviewed changes to each of the products listed as part of the VENTAK AV and VENTAK PRIZM DR/VR system and included this information in Section II when the change met the applicable criteria specified in 21 CFR 814.39 (b), for submission on a periodic reporting basis.

B. Recall and Safety Alert Update

1. On February 14, 2000, Guidant initiated a recall of a subset of the VENTAK PRIZM Models 1850/1851/1855/1856 devices. The action was initiated because of a random component failure in a subset of VENTAK PRIZM AICD devices. An assessment tool was distributed to the sales representatives and used to verify if a device is susceptible to the failure. FDA classified this as a "Recall" with reference number Z-468/470-0. A request for closure was submitted on April 25, 2000. The FDA closure letter is dated October 4, 2000.
2. On April 23, 2001, Guidant issued a letter to advise physicians of a situation that potentially affects a subset of VENTAK PRIZM 1 and VENTAK PRIZM HE devices, Models 1850, 1851, 1852, 1853, 1857 and 1858. Some of these devices have automatically switched to an integrated Safety Mode because of a rare interaction between the device and a specific memory component. No specific physician action was required. FDA classified this action as a "Recall" with reference number Z-613/618-1. A request for closure was submitted on June 25, 2001.

II. DEVICE CHANGES

This section describes those device changes that have been reported in accordance with 21 CFR 814.39(a), in a PMA supplemental submission, as well as those changes meeting the criteria specified in 21 CFR 814.39(b), for submission on a periodic reporting basis.

A. PMA Supplemental Submissions

The following supplemental changes have been made to the VENTAK PRIZM system (P960040) during the time period of this report.

- S17 Request approval for Version 4.0 to Application Software Model 2920 on PRx, AV and Application Software Model 2833 on AV II DDD. Submitted on July 24, 2000. Approved: December 13, 2000.
- S18 Request approval for ASTRID (Atrial Sensing to Reduce Inappropriate Defibrillation) claims associated with VENTAK AV. Submitted on September 15, 2000. Approval Pending.
- S19 Request approval (Real Time Review) for Torque Wrench Model 6942. Submitted on November 30, 2000. Approved: January 23, 2001.
- S20 Request approval (Real Time Review) for Design modification to the header on VENTAK PRIZM Model 1851, DR/HE Model 1853 and VR Models 1850/1852. Submitted on December 22, 2000. Approved: February 1, 2001.
- S21 Request approval for Version 3.3 to Application Software Model 2844 on VENTAK PRIZM. Submitted February 23, 2001. Approved: April 5, 2001.
- S22 Request approval (Real Time Review) for Diagnostic/Restoration Tool Version 1.5 for Memory Interaction Safety Mode Model 2725. Submitted May 23, 2001. Approved: June 13, 2001

B. Changes Not Included in a Submission

The following changes constitute minor manufacturing changes or other minor alterations to the device since the PMA filing date. These changes, which do not affect the safety and effectiveness of the device, do not require a PMA supplement per 21 CFR 814.39(a). These changes are therefore included in this post approval report as required in 21 CFR 814.39(b). Changes cover the period June 1, 2000 through May 31, 2001.

1. A manufacturing software test was modified to [REDACTED]

2. Guidant incorporated the capacitors [REDACTED]

[REDACTED] This is being done to standardize [REDACTED]

[REDACTED] All capacitors are 100% tested and meet all requirements of the capacitor specification. Device performance is unaffected by this change, and the device continues to meet the physical and functional performance requirements of Guidant device specification.

3. Guidant replaced the [REDACTED]

[REDACTED]
Device performance is unaffected by this change, and the device continues to meet the physical and functional performance requirements of Guidant device specification.

III. BIBLIOGRAPHY and SUMMARY of STUDIES

A. Unpublished Reports Regarding the Device

Guidant is aware of the following non-clinical investigations that involve pulse generator models in this report.

1. *Evaluation of Overdrive Dual-Chamber Pacing to Reduce the Incidence of Ventricular Tachyarrhythmias with the Guidant/CPI VENTAK AV AICD System.*

[REDACTED]

2. *Evaluation of Incremental Detection Enhancements and Delivery of Appropriate Therapy with the Guidant VENTAK AV System.*

[REDACTED]

3. *Comparison of Dual- and Single-Chamber Therapy Responsiveness Using the Guidant VENTAK® AICD Family.*

[REDACTED]

B. Reports in Scientific Literature Concerning the Device

Guidant conducted a literature search of English language articles in peer reviewed scientific/medical publications. Guidant searched files in the Medline database. Search descriptors were "implantable and defibrillator", biphasic cardioversion, and defibrillation shocks irrespective of manufacturer. The search did not specifically request accessories. The literature search includes articles dated between June 1, 2000 and May 31, 2001.

Reference List

1. Defibrillator/monitor/pacemakers. Health Devices. 2000 Sep; 29(9):301-34; ISSN: 0046-7022.
2. Keeping pace with cardiac implants. Johns Hopkins Med Lett Health After 50. 2000 Aug; 12(6):1-2; ISSN: 1042-1882.
3. Ahmad, M.; Bloomstein, L.; Roelke, M.; Bernstein, A. D., and Parsonnet, V. Patients' attitudes toward implanted defibrillator shocks. Pacing Clin Electrophysiol. 2000 Jun; 23(6):934-8; ISSN: 0147-8389.
4. Andraghetti, A. and Scalese, M. Safety and efficacy of low-energy cardioversion of 500 patients using two different techniques. Europace. 2001 Jan; 3(1):4-9; ISSN: 1099-5129.
5. Bansch, D.; Castrucci, M.; Bocker, D.; Breithardt, G., and Block, M. Ventricular tachycardias above the initially programmed tachycardia detection interval in patients with implantable cardioverter-defibrillators: incidence, prediction and significance. J Am Coll Cardiol. 2000 Aug; 36(2):557-65; ISSN: 0735-1097.
6. Blitzer, M. L.; Marieb, M. A., and Schoenfeld, M. H. Inability to communicate with ICDs: an underreported failure mode. Pacing Clin Electrophysiol. 2001 Jan; 24(1):13-5; ISSN: 0147-8389.
7. Brugada, J.; Brugada, R., and Brugada, P. Pharmacological and device approach to therapy of inherited cardiac diseases associated with cardiac arrhythmias and sudden death. J Electrocardiol. 2000; 33 Suppl41-7; ISSN: 0022-0736.
8. Brugada, P. Brugada syndrome: an electrocardiographic diagnosis not to be missed. Heart. 2000 Jul; 84(1):1-2; ISSN: 1468-201X.
9. Brugada, P.; Brugada, R., and Brugada, J. The Brugada syndrome. Curr Cardiol Rep. 2000 Nov; 2(6):507-14; ISSN: 1523-3782.
10. Butler, J. M. Brugada syndrome—the missed epidemic. J Accid Emerg Med. 2000 Nov; 17(6):425-8; ISSN: 1351-0622.

11. Capucci, A. Hospitalization in ICD recipients. *Eur Heart J.* 2000 Jul; 21(14):1123-4; ISSN: 0195-668X.
12. Chattipakorn, N.; Fotuhi, P. C., and Ideker, R. E. Pacing following shocks stronger than the defibrillation threshold: impact on defibrillation outcome. *J Cardiovasc Electrophysiol.* 2000 Sep; 11(9):1022-8; ISSN: 1045-3873.
13. Chua, J. D.; Wilkoff, B. L.; Lee, I.; Juratli, N.; Longworth, D. L., and Gordon, S. M. Diagnosis and management of infections involving implantable electrophysiologic cardiac devices. *Ann Intern Med.* 2000 Oct 17; 133(8):604-8; ISSN: 0003-4819.
14. Cohen, T. J.; Zadeh, H.; Fruauff, A., and Klein, J. Arrhythmogenic right ventricular dysplasia and the implantable cardioverter-defibrillator: a case report and review of the literature. *J Invasive Cardiol.* 2000 Aug; 12(8):422-4; ISSN: 1042-3931.
15. Connolly, S. J.; Krahn, A., and Klein, G. Long term management of the survivor of ventricular fibrillation or sustained ventricular tachycardia. *Can J. Cardiol.* 2000 Jun; 16 Suppl C20C-2C; ISSN: 0828-282X.
16. Costeas, X. F.; Link, M. S.; Foote, C. B.; Homoud, M. K.; Wang, P. J., and Estes, N. A. Predictors of ventricular tachycardia recurrence in 100 patients receiving tiered therapy defibrillators. *Clin Cardiol.* 2000 Nov; 23(11):852-6; ISSN: 0160-9289.
17. Costeas, X. F.; Strembelas, P. G.; Markou, D. X.; Stefanadis, C. I., and Toutouzas, P. K. Subpectoral cardioverter-defibrillator implantation using a lateral approach. *J Interv Card Electrophysiol.* 2000 Dec; 4(4):611-9; ISSN: 1383-875X.
18. De Sutter, J.; Tavernier, R.; Van de Wiele, C.; Kazmierczak, J.; De Buyzere, M.; Jordaens, L.; Clement, D. L., and Dierckx, R. A. Infarct size and recurrence of ventricular arrhythmias after defibrillator implantation. *Eur J Nucl Med.* 2000 Jul; 27(7):807-15; ISSN: 0340-6997.
19. Delacretaz, E.; Schlaepfer, J.; Metzger, J.; Fromer, M., and Kappenberger, L. Evidence rather than costs must guide use of the implantable cardioverter defibrillator. *Am J Cardiol.* 2000 Nov 2; 86(9 Suppl 1):K52-K57; ISSN: 0002-9149.
20. Dickens, B. M. Chapter 5. Legal implications of ICD therapy. *Can J Cardiol.* 2000 Oct; 16(10):1319-24; ISSN: 0828-282X.
21. Dickerson, S. S.; Flaig, D. M., and Kennedy, M. C. Therapeutic connection: help seeking on the Internet for persons with implantable cardioverter defibrillators. *Heart Lung.* 2000 Jul-2000 Aug 31; 29(4):248-55; ISSN: 0147-9563.

22. ---. Dual chamber arrhythmia detection in the implantable cardioverter defibrillator. *J Cardiovasc Electrophysiol*. 2000 Oct; 11(10):1105-15; ISSN: 1045-3873.
23. ---. Importance of the atrial channel for ventricular arrhythmia therapy in the dual chamber implantable cardioverter defibrillator. *J Cardiovasc Electrophysiol*. 2000 Dec; 11(12):1309-19; ISSN: 1045-3873.
24. Doherty, J. U.; Fuchs, S., and Tecce, M. A. Ventricular arrhythmias. Preventing sudden death with drugs and ICD devices. *Geriatrics*. 2000 Aug; 55(8):26-8, 31-2, 35-6, ISSN:0016-867X.
25. ---. Combination ICD and drug treatments-best options. *Resuscitation*. 2000 Aug 1; 45(3):S3-6; ISSN: 0300-9572.
26. Duru, F.; Bauersfeld, U., and Candinas, R. Avoiding inappropriate ventricular tachycardia detection due to T-wave oversensing in an implantable cardioverter defibrillator: a novel application of the electrogram width criterion. *Europace*. 2001 Jan; 3(1):80-4; ISSN: 1099-5129.
27. Duru, F.; Bauersfeld, U.; Rahn-Schonbeck, M., and Candinas, R. Morphology discriminator feature for enhanced ventricular tachycardia discrimination in implantable cardioverter defibrillators. *Pacing Clin Electrophysiol*. 2000 Sep; 23(9):1365-74; ISSN: 0147-8389.
28. Duru, F.; Buchi, S.; Klaghofer, R.; Mattmann, H.; Sensky, T.; Buddeberg, C., and Candinas, R. How different from pacemaker patients are recipients of implantable cardioverter-defibrillators with respect to psychosocial adaptation, affective disorders, and quality of life? *Heart*. 2001 Apr; 85(4):375-9; ISSN: 1468-201X.
29. Eckart, R. E.; Kinney, K. G.; Belnap, C. M., and Le, T. D. Ventricular fibrillation refractory to automatic internal cardiac defibrillator in Fabry's disease. Review of cardiovascular manifestations. *Cardiology*. 2000; 94(3):208-12; ISSN: 0008-6312.
30. El-Gamal, H. M.; Dufresne, R. G., and Saddler, K. Electrosurgery, pacemakers and ICDs: a survey of precautions and complications experienced by cutaneous surgeons. *Dermatol Surg*. 2001 Apr; 27(4):385-90; ISSN: 1076-0512.
31. Ellenbogen, K. A.; Edel, T.; Moore, S.; Higgins, S.; Pacifico, A.; Wilber, D.; Wood, M. A.; Rogers, R.; Dahn, A., and Zhu, A. A prospective randomized-controlled trial of ventricular fibrillation detection time in a DDDR ventricular defibrillator. Ventak AV II DR Study Investigators. *Pacing Clin Electrophysiol*. 2000 Aug; 23(8):1268-72; ISSN: 0147-8389.
32. Exner, D. V.; Sheldon, R. S.; Pinski, S. L.; Kron, J., and Hallstrom, A. Do baseline characteristics accurately discriminate between patients likely versus unlikely to benefit from implantable defibrillator therapy? Evaluation of the

12

Canadian implantable defibrillator study implantable cardioverter defibrillatory efficacy score in the antiarrhythmics versus implantable defibrillators trial. *Am Heart J*. 2001 Jan; 141(1):99-104; ISSN: 0002-8703.

33. Fridlund, B.; Lindgren, E. C.; Ivarsson, A.; Jinhage, B. M.; Boise, K.; Flemme, I.; Sandstedt, B., and Martensson, J. Patients with implantable cardioverter-defibrillators and their conceptions of the life situation: a qualitative analysis. *J Clin Nurs*. 2000 Jan; 9(1):37-45; ISSN: 0962-1067.
34. Friedman, P. A.; Glikson, M., and Stanton, M. S. Defibrillator challenges for the new millennium: the marriage of device and patient-making and maintaining a good match. *J Cardiovasc Electrophysiol*. 2000 Jun; 11(6):697-709; ISSN: 1045-3873.
35. Gatzoulis, K.; Protonotarios, N.; Anastasakis, A.; Tsatsopoulou, A.; Vlasseros, J.; Gialafos, J., and Toutouzas, P. Implantable defibrillator therapy in Naxos disease. *Pacing Clin Electrophysiol*. 2000 Jul; 23(7):1176-8; ISSN: 0147-8389.
36. Geist, M.; Newman, D.; Greene, M.; Paquette, M., and Dorian, P. Permanent explanation of implantable cardioverter defibrillators. *Pacing Clin Electrophysiol*. 2000 Dec; 23(12):2024-9; ISSN: 0147-8389.
37. Glikson, M. and Friedman, P. A. The implantable cardioverter defibrillator. *Lancet*. 2001 Apr 7; 357(9262):1107-17; ISSN: 0140-6736.
38. Godemann, F.; Ahrens, B.; Behrens, S.; Berthold, R.; Gandor, C.; Lampe, F., and Linden, M. Classic conditioning and dysfunctional cognitions in patients with panic disorder and agoraphobia treated with an implantable cardioverter/defibrillator. *Psychosom Med*. 2001 Mar-2001 Apr 30; 63(2):231-8; ISSN: 0033-3174.
39. Gold, M. R. ICD therapy in the new millennium. *Cardiol Clin*. 2000 May; 18(2):375-89; ISSN: 0733-8651.
40. Gold, M. R. and Nisam, S. Primary prevention of sudden cardiac death with implantable cardioverter defibrillators: lessons learned from MADIT and MUSTT. *Pacing Clin Electrophysiol*. 2000 Nov; 23(11 Pt 2):1981-5; ISSN: 0147-8389.
Notes: CORPORATE NAME: Multicenter Automatic Defibrillator Implantation Trial. Multicenter Unsustained Tachycardia Trial.
41. Goldschlager, N.; Epstein, A.; Friedman, P.; Gang, E.; Krol, R., and Olshansky, B. Environmental and drug effects on patients with pacemakers and implantable cardioverter/defibrillators: a practical guide to patient treatment. *Arch Intern Med*. 2001 Mar 12; 161(5):649-55; ISSN: 0003-9926.
Notes: CORPORATE NAME: North American Society of Pacing and Electrophysiology (NASPE) Practice Guideline Committee.

42. Gollob, M. H. and Seger, J. J. Current status of the implantable cardioverter-defibrillator. *Chest*. 2001 Apr; 119(4):1210-21; ISSN: 0012-3692.
43. Gurevitz, O.; Yaacoby, E.; Segal, E.; Perel, A.; Eldar, M.; Kaplinsky, E., and Glikson, M. Effect of implantable cardioverter-defibrillator shocks on QT dispersion. *Am J Cardiol*. 2000 Nov 15; 86(10):1146-8, A9; ISSN: 0002-9149.
44. Hallstrom, A. P.; Anderson, J. L.; Cobb, L. A.; Friedman, P. L.; Herre, J. M.; Klein, R. C.; McAnulty, J., and Steinberg, J. S. Advantages and disadvantages of trial designs: a review of analysis methods for ICD studies. AVID Investigators. *Pacing Clin Electrophysiol*. 2000 Jun; 23(6):1029-38; ISSN: 0147-8389.
45. Hallstrom, A. P.; McAnulty, J. H.; Wilkoff, B. L.; Follmann, D.; Raitt, M. H.; Carlson, M. D.; Gillis, A. M.; Shih, H. T.; Powell, J. L.; Duff, H., and Halperin, B. D. Patients at lower risk of arrhythmia recurrence: a subgroup in whom implantable defibrillators may not offer benefit. Antiarrhythmics Versus Implantable Defibrillator (AVID) Trial Investigators. *J Am Coll Cardiol*. 2001 Mar 15; 37(4):1093-9; ISSN: 0735-1097.
Notes: CORPORATE NAME: Antiarrhythmics Versus Implantable Defibrillator (AVID) Trial Investigators.
46. Hauser, R.; Hayes, D.; Parsonnet, V.; Furman, S.; Epstein, A.; Hayes, J.; Saksena, S.; Irwin, M.; Almquist, A.; Cannom, D.; Gross, J., and Kallinen, L. Feasibility and initial results of an Internet-based pacemaker and ICD pulse generator and lead registry. *Pacing Clin Electrophysiol*. 2001 Jan; 24(1):82-7; ISSN: 0147-8389.
47. Hazekamp, M. G.; Blom, N. A.; Schoof, P. H.; Schalij, M. J., and Dion, R. A. Implantation of cardioverter device in young children: the perirenal approach. *Ann Thorac Surg*. 2001 Apr; 71(4):1382-3; ISSN: 0003-4975.
48. Hegbom, F.; Hoff, P. I.; Oie, B.; Folling, M.; Zeijlemaker, V.; Lindemans, F., and Ohm, O. J. RV function in stable and unstable VT: is there a need for hemodynamic monitoring in future defibrillators? *Pacing Clin Electrophysiol*. 2001 Feb; 24(2):172-82; ISSN: 0147-8389.
49. Hintringer, F. Highest number of ICD implantations per million inhabitants from 1993-1998 in Austria. *Pacing Clin Electrophysiol*. 2001 Jan; 24(1):129. ISSN: 0147-8389.
50. Hoffmaster, B. Chapter 4. The ethics of setting limits on ICD therapy. *Can J Cardiol*. 2000 Oct; 16(10):1313-8; ISSN: 0828-282X.
51. Hohnloser, S. H.; Connolly, S. J.; Kuck, K. H.; Dorian, P.; Fain, E.; Hampton, J. R.; Hatala, R.; Pauly, A. C.; Roberts, R. S.; Themeles, E., and Gent, M. The defibrillator in acute myocardial infarction trial (DINAMIT): study protocol.

Am Heart J. 2000 Nov; 140(5):735-9; ISSN: 0002-8703.

52. Ideker, R. E.; Chattipakorn, T. N., and Gray, R. A. Defibrillation mechanisms: the parable of the blind men and the elephant. *J Cardiovasc Electrophysiol*. 2000 Sep; 11(9):1008-13; ISSN: 1045-3873.
53. Kaczmarek, R. G.; Beaulieu, M. D., and Kessler, L. G. Medical device tracking: results of a case study of the implantable cardioverter defibrillator. *Am J Cardiol*. 2000 Mar 1; 85(5):588-92; ISSN: 0002-9149.
54. Kakishita, M.; Kurita, T.; Matsuo, K.; Taguchi, A.; Suyama, K.; Shimizu, W.; Aihara, N.; Kamakura, S.; Yamamoto, F.; Kobayashi, J.; Kosakai, Y., and Ohe, T. Mode of onset of ventricular fibrillation in patients with Brugada syndrome detected by implantable cardioverter defibrillator therapy. *J Am Coll Cardiol*. 2000 Nov 1; 36(5):1646-53; ISSN: 0735-1097.
55. Karagueuzian, H. S. and Chen, P. S. Fibrillation and defibrillation: the odd couple? *J Cardiovasc Electrophysiol*. 2000 Jun; 11(6):642-4; ISSN: 1045-3873.
56. Kesek, M.; Saraj, A., and Blomstrom-Lundqvist, C. The value of predischage ICD tests in patients with a successful peroperative test. *Scand Cardiovasc J*. 2000; 34(2):130-3; ISSN: 1401-7431.
57. Kolb, C.; Zrenner, B., and Schmitt, C. Incidence of electromagnetic interference in implantable cardioverter defibrillators. *Pacing Clin Electrophysiol*. 2001 Apr; 24(4 Pt 1):465-8; ISSN: 0147-8389.
58. Korte, T.; Jung, W.; Ostermann, G.; Wolpert, C.; Spehl, S.; Esmailzadeh, B., and Luderitz, B. Hospital readmission after transvenous cardioverter/defibrillator implantation; a single centre study. *Eur Heart J*. 2000 Jul; 21(14):1186-91; ISSN: 0195-668X.
59. Krishnan, S. C. and Josephson, M. E. Surgery for postinfarction ventricular tachycardia: is it obsolete? *Pacing Clin Electrophysiol*. 2000 Aug; 23(8):1295-301; ISSN: 0147-8389.
60. Kron, J.; Herre, J.; Renfro, E. G.; Rizo-Patron, C.; Raitt, M.; Halperin, B.; Gold, M.; Goldner, B.; Wathen, M.; Wilkoff, B.; Olarte, A., and Yao, Q. Lead- and device-related complications in the antiarrhythmics versus implantable defibrillators trial. *Am Heart J*. 2001 Jan; 141(1):92-8; ISSN: 0002-8703.
61. Kuck, K. H.; Cappato, R.; Siebels, J., and Ruppel, R. Randomized comparison of antiarrhythmic drug therapy with implantable defibrillators in patients resuscitated from cardiac arrest : the Cardiac Arrest Study Hamburg (CASH). *Circulation*. 2000 Aug 15; 102(7):748-54; ISSN: 1524-4539.
62. Lewin, R. J.; Frizelle, D. J., and Kaye, G. C. A rehabilitative approach to patients with internal cardioverter-defibrillators. *Heart*. 2001 Apr; 85(4):371-2;

ISSN: 1468-201X.

63. Lin, Y. C.; Chang, W. T.; Chen, W. J.; Lai, L. P.; Lin, J. L., and Wu, M. H. Control of idiopathic ventricular fibrillation by implantable cardioverter-defibrillator in a child who survived sudden death. *J Formos Med Assoc.* 2000 Jul; 99(7):576-9; ISSN: 0929-6646.
64. Linde, C. Quality-of-life in pacemaker and implantable cardioverter defibrillator recipients. *Pacing Clin Electrophysiol.* 2000 Jun; 23(6):931-3; ISSN: 0147-8389.
65. Lombardi, F.; Porta, A.; Marzegalli, M.; Favale, S.; Santini, M.; Vincenti, A., and De Rosa, A. Heart rate variability patterns before ventricular tachycardia onset in patients with an implantable cardioverter defibrillator. Participating Investigators of ICD-HRV Italian Study Group. *Am J Cardiol.* 2000 Nov 1; 86(9):959-63; ISSN: 0002-9149.
Notes: CORPORATE NAME: Implantable Cardioverter Defibrillator-Heart
66. Lubinski, A.; Lewicka-Nowak, E.; Krolak, T.; Kempa, M.; Bielawska, B.; Wilczek, R., and Swiatecka, G. Implantation and follow-up of ICD leads implanted in the right ventricular outflow tract. *Pacing Clin Electrophysiol.* 2000 Nov; 23(11 Pt 2):1996-8; ISSN: 0147-8389.
67. Luceri, R. M. Initial clinical experience with a dual chamber rate responsive implantable cardioverter defibrillator. *Pacing Clin Electrophysiol.* 2000 Nov; 23(11 Pt 2):1986-8; ISSN: 0147-8389.
Notes: CORPORATE NAME: Photon DR Clinical Investigators.
68. Mann, D. E.; Otto, L.; Kelly, P. A., and Reiter, M. J. Effect of sensing system on the incidence of myopotential oversensing during bradycardia pacing in implantable cardioverter-defibrillators. *Am J Cardiol.* 2000 Jun 1; 85(11):1380-2; ISSN: 0002-9149.
69. Manolis, A. S.; Chiladakis, J.; Maounis, T. N.; Vassilikos, V., and Cokkinos, D. V. Two-coil versus single-coil transvenous cardioverter defibrillator systems: comparative data. *Pacing Clin Electrophysiol.* 2000 Nov; 23(11 Pt 2):1999-2002; ISSN: 0147-8389.
70. Mazur, A.; Anderson, M. E.; Bonney, S., and Roden, D. M. Pause-dependent polymorphic ventricular tachycardia during long-term treatment with dofetilide: a placebo-controlled, implantable cardioverter-defibrillator-based evaluation. *J Am Coll Cardiol.* 2001 Mar 15; 37(4):1100-5; ISSN: 0735-1097.
71. Mazur, A.; Wang, L.; Anderson, M. E.; Yee, R.; Theres, H.; Pearson, A.; Olson, W., and Wathen, M. Functional similarity between electrograms recorded from an implantable cardioverter defibrillator emulator and the surface electrocardiogram. *Pacing Clin Electrophysiol.* 2001 Jan; 24(1):34-40;

16

ISSN: 0147-8389.

72. Michaud, G. F.; Pelosi, F. Jr; Noble, M. D.; Knight, B. P.; Morady, F., and Strickberger, S. A. A randomized trial comparing heparin initiation 6 h or 24 h after pacemaker or defibrillator implantation. *J Am Coll Cardiol.* 2000 Jun; 35(7):1915-8; ISSN: 0735-1097.
73. Mitchell, L. B. and Roy, D. Treatment approaches for patients with sustained VT/VF. *Can J Cardiol.* 2000 Jun; 16 Suppl C23C-5C; ISSN: 0828-282X.
74. Morris, M. M.; KenKnight, B. H., and Lang, D. J. Detection of atrial arrhythmia for cardiac rhythm management by implantable devices. *J Electrocardiol.* 2000; 33 Suppl133-9; ISSN: 0022-0736.
75. Mouchawar, G.; Kroll, M.; Val-Mejias, J. E.; Schwartzman, D.; McKenzie, J.; Fitzgerald, D.; Prater, S.; Katcher, M.; Fain, E., and Syed, Z. ICD waveform optimization: a randomized, prospective, pair-sampled multicenter study. *Pacing Clin Electrophysiol.* 2000 Nov; 23(11 Pt 2):1992-5; ISSN: 0147-8389.
76. Myerburg, R. J. Life-threatening arrhythmias and therapy. *Resuscitation.* 2000 Aug 1; 45(3):S1-2; ISSN: 0300-9572.
77. Naik, A. M. and Peter, C. T. Are implantable cardioverter- defibrillators cost-effective? *Curr Cardiol Rep.* 2000 Jul; 2(4):341-5; ISSN: 1523-3782.
78. Nair, L. A.; Ostrow, E. L., and Beck, B. ICD hardware failure associated with multiple internal shocks. *Pacing Clin Electrophysiol.* 2000 May; 23(5):910-3; ISSN: 0147-8389.
79. Neel, M. Posttraumatic stress symptomatology in patients with automatic implantable cardioverter defibrillators: nature and intervention. *Int J Emerg Ment Health.* 2000 Fall; 2(4):259-63; ISSN: 1522-4821.
80. Noorani, H. Z. Chapter 2. Guidelines and recommendations for ICD therapy. *Can J Cardiol.* 2000 Oct; 16(10):1303-6; ISSN: 0828-282X.
81. Nordberg, M. Cardiac care 2000. A shocking discovery. *Emerg Med Serv.* 2000 Sep; 29(9):67, 69-70, 72 passim; ISSN: 0094-6575.
82. O'Brien, B. J.; Connolly, S. J.; Goeree, R.; Blackhouse, G.; Willan, A.; Yee, R.; Roberts, R. S., and Gent, M. Cost-effectiveness of the implantable cardioverter-defibrillator: results from the Canadian Implantable Defibrillator Study (CIDS). *Circulation.* 2001 Mar 13; 103(10):1416-21; ISSN: 1524-4539.
83. O'Callaghan, P. A. Prophylactic defibrillator implantation in patients with coronary artery disease. *Curr Cardiol Rep.* 2000 Jul; 2(4):311-3; ISSN: 1523-3782.

84. Obwegeser, A. A.; Uitti, R. J.; Turk, M. F.; Wszolek, U. M.; Flipse, T. R.; Smallridge, R. C.; Witte, R. J., and Wharen, R. E. Jr. Simultaneous thalamic deep brain stimulation and implantable cardioverter-defibrillator. *Mayo Clin Proc.* 2001 Jan; 76(1):87-9; ISSN: 0025-6196.
85. Ocampo, C. M. Living with an implantable cardioverter defibrillator: impact on the patient, family, and society. *Nurs Clin North Am.* 2000 Dec; 35(4):1019-30; ISSN: 0029-6465.
86. Page, R. L. Effects of antiarrhythmic medication on implantable cardioverter-defibrillator function. *Am J Cardiol.* 2000 Jun 15; 85(12):1481-5; ISSN: 0002-9149.
87. Parkes, J.; Bryant, J., and Milne, R. Implantable cardioverter defibrillators: arrhythmias. A rapid and systematic review. *Health Technol Assess.* 2000; 4(26):1-69; ISSN: 1366-5278.
88. Passman, R. and Kadish, A. Polymorphic ventricular tachycardia, long Q-T syndrome, and torsades de pointes. *Med Clin North Am.* 2001 Mar; 85(2):321-41; ISSN: 0025-7125.
89. Pavia, S. and Wilkoff, B. L. Preventing sudden death in coronary cardiomyopathy: implantable defibrillators lead the way. *Cleve Clin J Med.* 2001 Feb; 68(2):113, 118, 120, 121, 125, 129; ISSN: 0891-1150.
90. Peters, R. W. and Gold, M. R. Implantable cardiac defibrillators. *Med Clin North Am.* 2001 Mar; 85(2):343-67, xi; ISSN: 0025-7125.
91. Peters, R. W.; Zhang, X., and Gold, M. R. Clinical predictors and efficacy of antitachycardia pacing in patients with implantable cardioverter defibrillators: the importance of the patient's sex. *Pacing Clin Electrophysiol.* 2001 Jan; 24(1):70-4; ISSN: 0147-8389.
92. Pfeiffer, D.; Mende, M., and Hagendorff, A. Implantable dual-chamber cardioverter-defibrillator-pacemaker. *Curr Cardiol Rep.* 2000 Jul; 2(4):335-40; ISSN: 1523-3782.
93. Pinski, S. L. Emergencies related to implantable cardioverter-defibrillators. *Crit Care Med.* 2000 Oct; 28(10 Suppl):N174-80; ISSN: 0090-3493.
94. ---. Inappropriate pacing due to autoperpetuation of the ventricular rate stabilization algorithm: a manifestation of T wave oversensing by ICDs. *Pacing Clin Electrophysiol.* 2000 Sep; 23(9):1446-7; ISSN: 0147-8389.
95. Pinski, S. L. and Bredikis, A. J. Defibrillator implantation via the iliac vein. *Pacing Clin Electrophysiol.* 2000 Aug; 23(8):1315-7; ISSN: 0147-8389.
96. Pinski, S. L. and Trohman, R. G. Permanent pacing via implantable defibrillators.

18

- Pacing Clin Electrophysiol. 2000 Nov; 23(11 Pt 1):1667-82; ISSN: 0147-8389.
97. Prystowsky, E. ICD replacement: a family-based medicine approach. Pacing Clin Electrophysiol. 2000 Dec; 23(12):2022-3; ISSN: 0147-8389.
98. Prystowsky, E. N. and Nisam, S. Prophylactic implantable cardioverter defibrillator trials: MUSTT, MADIT, and beyond. Multicenter Unsustained Tachycardia Trial. Multicenter Automatic Defibrillator Implantation Trial. Am J Cardiol. 2000 Dec 1; 86(11):1214-5, A5; ISSN: 0002-9149.
99. Raj, S. R.; Mitchell, L. B., and Sheldon, R. S. Implantable cardioverter-defibrillators in hypertrophic cardiomyopathy. N Engl J Med. 2000 Jul 6; 343(1):68-9. ISSN: 0028-4793.
100. Saeed, M.; Link, M. S.; Mahapatra, S.; Mouded, M.; Tzeng, D.; Jung, V.; Contreras, R.; Swygman, C.; Homoud, M.; Estes, N. A. 3rd, and Wang, P. J. Analysis of intracardiac electrograms showing monomorphic ventricular tachycardia in patients with implantable cardioverter-defibrillators. Am J Cardiol. 2000 Mar 1; 85(5):580-7; ISSN: 0002-9149.
101. Saksena, S. Implantable defibrillators in the third millennium: increasingly relegated to a standby role? J Am Coll Cardiol. 2000 Sep; 36(3):828-31; ISSN: 0735-1097.
102. Schaer, B. and Osswald, S. Methods of minimizing inappropriate implantable cardioverter-defibrillator shocks. Curr Cardiol Rep. 2000 Jul; 2(4):346-52; ISSN: 1523-3782.
103. Schluter, T.; Baum, H.; Plewan, A., and Neumeier, D. Effects of implantable cardioverter defibrillator implantation and shock application on biochemical markers of myocardial damage. Clin Chem. 2001; 47(3):459-63; ISSN: 0009-9147.
104. Schmitt, H.; Hurst, T.; Coch, M.; Killat, H.; Wunn, B., and Waldecker, B. Nonsustained, asymptomatic ventricular tachycardia in patients with coronary artery disease: prognosis and incidence of sudden death of patients who are noninducible by electrophysiological testing. Pacing Clin Electrophysiol. 2000 Aug; 23(8):1220-5; ISSN: 0147-8389.
105. Schoenfeld, M. H. Deciding against defibrillator replacement: second-guessing the past? Pacing Clin Electrophysiol. 2000 Dec; 23(12):2019-21; ISSN: 0147-8389.
106. Schreiber, C.; Mehmanesch, H.; Kolb, C.; Schmitt, C., and Lange, R. Modified implantation of a transvenous defibrillator in a patient after tricuspid valve replacement. Pacing Clin Electrophysiol. 2000 Nov; 23(11 Pt 1):1698-9; ISSN: 0147-8389.

107. Schulte, B.; Sperzel, J.; Schwarz, T.; Pitschner, H. F.; Strupp, G., and Neuzner, J. Detection of ventricular fibrillation in implantable defibrillators with automatic gain control amplifiers: effects of programming sensitivity. *Europace*. 2000 Apr; 2(2):160-2; ISSN: 1099-5129.
108. Sears, S. F.; Todaro, J. F.; Urizar, G.; Lewis, T. S.; Sirois, B.; Wallace, R.; Sotile, W.; Curtis, A. B., and Conti, J. B. Assessing the psychosocial impact of the ICD: a national survey of implantable cardioverter defibrillator health care providers. *Pacing Clin Electrophysiol*. 2000 Jun; 23(6):939-45; ISSN: 0147-8389.
109. Shellock, F. G. MR imaging and electronically activated devices. *Radiology*. 2001 Apr; 219(1):294-5. ISSN: 0033-8419.
110. Shivkumar, K.; Feliciano, Z.; Boyle, N. G., and Wiener, I. Intradevice interaction in a dual chamber implantable cardioverter defibrillator preventing ventricular tachyarrhythmia detection. *J Cardiovasc Electrophysiol*. 2000 Nov; 11(11):1285-8; ISSN: 1045-3873.
111. Simpson, C. S.; Klein, G. J., and Hoffmaste, B. Expensive medical technologies and "indication extrapolation": the case of implantable cardioverter-defibrillators. *Am Heart J*. 2000 Sep; 140(3):419-22; ISSN: 0002-8703.
112. Singh, S. N.; Karasik, P.; Hafley, G. E.; Pieper, K. S.; Lee, K. L.; Wyse, D. G., and Buxton, A. E. Electrophysiologic and clinical effects of angiotensin-converting enzyme inhibitors in patients with prior myocardial infarction, nonsustained ventricular tachycardia, and depressed left ventricular function. MUSTT Investigators. Multicenter UnSustained Tachycardia Trial. *Am J Cardiol*. 2001 Mar 15; 87(6):716-20; ISSN: 0002-9149. Notes: CORPORATE NAME: MUSTT Investigators. Multicenter UnSustained Tachycardia Trial.
113. Smith, E. R. Implantable cardioverter/defibrillator therapy: the meaning of terms. *Can J Cardiol*. 2000 Oct; 16(10):1325-6; ISSN: 0828-282X.
114. SoRelle, R. Defibrillator disagreement. *Circulation*. 2000 Oct 3; 102(14):E9028-9; ISSN: 1524-4539.
115. Splett, V.; Trusty, J. M.; Hayes, D. L., and Friedman, P. A. Determination of pacing capture in implantable defibrillators: benefit of evoked response detection using RV coil to can vector. *Pacing Clin Electrophysiol*. 2000 Nov; 23(11 Pt 1):1645-50; ISSN: 0147-8389.
116. Sra, J. and Akhtar, M. Inappropriate shock delivery by implantable defibrillators with dual chamber pacing during nonsustained ventricular tachycardia in patients with heart block. *Pacing Clin Electrophysiol*. 2000 Jun; 23(6):1054-6; ISSN: 0147-8389.

117. Subramanian, V. A. and Patel, N. U. Canadian implantable defibrillator study: what does it mean after the antiarrhythmics versus implantable defibrillators trial? *Curr Cardiol Rep.* 2000 Nov; 2(6):489-91; ISSN: 1523-3782.
118. Sun, S.; Klouche, K.; Tang, W., and Weil, M. H. The effects of biphasic and conventional monophasic defibrillation on postresuscitation myocardial function. *J Am Coll Cardiol.* 2001 May; 37(6):1753-4. ISSN: 0735-1097.
119. Sweeney, M. O.; Natale, A.; Volosin, K. J.; Swerdlow, C. D.; Baker, J. H., and Degroot, P. Prospective randomized comparison of 50%/50% versus 65%/65% tilt biphasic waveform on defibrillation in humans. *Pacing Clin Electrophysiol.* 2001 Jan; 24(1):60-5; ISSN: 0147-8389.
120. Talajic, M.; Mitchell, L. B., and Hadjis, T. The role of prophylactic implantable cardioverter defibrillators to prevent sudden cardiac death. *Can J Cardiol.* 2000 Jun; 16 Suppl C26C-8C; ISSN: 0828-282X.
121. Tavernier, R.; Fonteyne, W.; Vandewalle, V.; de Sutter, J., and Gevaert, S. Use of an implantable cardioverter defibrillator in a patient with two implanted neurostimulators for severe Parkinson's disease. *Pacing Clin Electrophysiol.* 2000 Jun; 23(6):1057-9; ISSN: 0147-8389.
122. Tavernier, R.; Gevaert, S.; De Sutter, J.; De Clercq, A.; Rottiers, H.; Jordaens, L., and Fonteyne, W. Long term results of cardioverter-defibrillator implantation in patients with right ventricular dysplasia and malignant ventricular tachyarrhythmias. *Heart.* 2001 Jan; 85(1):53-6; ISSN: 1468-201X.
123. Taylor, E.; Berger, R.; Hummel, J. D.; Dinerman, J. L.; Kenknight, B.; Arria, A. M.; Tomaselli, G., and Calkins, H. Analysis of the pattern of initiation of sustained ventricular arrhythmias in patients with implantable defibrillators. *J Cardiovasc Electrophysiol.* 2000 Jul; 11(7):719-26; ISSN: 1045-3873.
124. Theuns, D.; Kimman, G. J., and Jordaens, L. Images in cardiology: Apparent induction of ventricular tachycardia after "appropriate pacing" by an implantable dual chamber defibrillator: confusing ICD electrograms. *Heart.* 2000 Jul; 84(1):36; ISSN: 1468-201X.
125. Thomas, S. A.; Friedmann, E., and Kelley, F. J. Living with an implantable cardioverter-defibrillator: a review of the current literature related to psychosocial factors. *AACN Clin Issues.* 2001 Feb; 12(1):156-63; ISSN: 1079-0713.
126. Tovar, O. H. and Jones, J. L. Electrophysiologic deterioration after one-minute fibrillation increases relative biphasic defibrillation efficacy. *J Cardiovasc Electrophysiol.* 2000 Jun; 11(6):645-51; ISSN: 1045-3873.
127. Tresch, D. D. Evaluation and management of cardiac arrhythmias in the elderly.

el

Med Clin North Am. 2001 Mar; 85(2):527-50, xii; ISSN: 0025-7125.

128. Trohman, R. G. and Parrillo, J. E. Direct current cardioversion: indications, techniques, and recent advances. Crit Care Med. 2000 Oct; 28(10 Suppl):N170-3; ISSN: 0090-3493.
129. Unterberg, C.; Stevens, J.; Vollmann, D.; Hasenfuss, G., and Buchwald, A. B. Long-term clinical experience with the EGM width detection criterion for differentiation of supraventricular and ventricular tachycardia in patients with implantable cardioverter defibrillators. Pacing Clin Electrophysiol. 2000 Nov; 23(11 Pt 1):1611-7; ISSN: 0147-8389.
130. Vanhees, L.; Schepers, D.; Heidbuchel, H.; Defoor, J., and Fagard, R. Exercise performance and training in patients with implantable cardioverter-defibrillators and coronary heart disease. Am J Cardiol. 2001 Mar 15; 87(6):712-5; ISSN: 0002-9149.
131. Viskin, S. and Fish, R. Prevention of ventricular arrhythmias in the congenital long QT syndrome. Curr Cardiol Rep. 2000 Nov; 2(6):492-7; ISSN: 1523-3782.
132. Weber, M.; Block, M.; Bansch, D.; Castrucci, M.; Gradaus, R.; Schriever, C.; Breithardt, G., and Bocker, D. Antitachycardia pacing for rapid VT during ICD charging: a method to prevent ICD shocks. Pacing Clin Electrophysiol. 2001 Mar; 24(3):345-51; ISSN: 0147-8389.
133. Weigner, M. J. and Buxton, A. E. Nonsustained ventricular tachycardia. A guide to the clinical significance and management. Med Clin North Am. 2001 Mar; 85(2):305-20, x; ISSN: 0025-7125.
134. Wessel, N.; Ziehmman, C.; Kurths, J.; Meyerfeldt, U.; Schirdewan, A., and Voss, A. Short-term forecasting of life-threatening cardiac arrhythmias based on symbolic dynamics and finite-time growth rates. Phys Rev E Stat Phys Plasmas Fluids Relat Interdiscip Topics. 2000 Jan; 61(1):733-9; ISSN: 1063-651X.
135. Wickline, S. A. and Fischer, K. C. Can infections be imaged in implanted devices? ASAIO J. 2000 Nov-2000 Dec 31; 46(6):S80-1; ISSN: 1058-2916.
136. Wilkoff, B. L.; Kuhlkamp, V.; Volosin, K.; Ellenbogen, K.; Waldecker, B.; Kacet, S.; Gillberg, J. M., and DeSouza, C. M. Critical analysis of dual-chamber implantable cardioverter-defibrillator arrhythmia detection : results and technical considerations. Circulation. 2001 Jan 23; 103(3):381-6; ISSN: 1524-4539.
137. Wood, M. A. and Ellenbogen, K. A. Initiation of spontaneous ventricular tachycardia: from spark to fire. J Cardiovasc Electrophysiol. 2000 Jul; 11(7):727-9; ISSN: 1045-3873.

138. Yamamoto, T.; Katayama, Y.; Fukaya, C.; Kurihara, J.; Oshima, H., and Kasai, M. Thalamotomy caused by cardioversion in a patient treated with deep brain stimulation. *Stereotact Funct Neurosurg.* 2000; 74(2):73-82; ISSN: 1011-6125.
139. Yamanouchi, Y.; Fishler, M. G.; Mowrey, K. A.; Wilkoff, B. L.; Mazgalev, T. N., and Tchou, P. J. New approach to biphasic waveforms for internal defibrillation: fully discharging capacitors. *J Cardiovasc Electrophysiol.* 2000 Aug; 11(8):907-12; ISSN: 1045-3873.

IV. SYSTEM EXPERIENCE

The VENTAK AV/VR PRIZM System Experience Report is submitted in response to requirements specified in the Conditions of Approval and is found in Appendix B of this report.

V. LIFE TESTING RESULTS

In accordance with the conditions of approval for PMA 960040/S9, the capacitor life testing report for VENTAK VR is found in Appendix C of this report.

APPENDIX A: CONDITIONS OF APPROVAL

Issued: 3-4-98

CONDITIONS OF APPROVAL
FOR IMPLANTABLE DEFIBRILLATORS AND PROGRAMMERS

APPROVED LABELING. As soon as possible, and before commercial distribution of your device, submit three copies of an amendment to this PMA submission with copies of all approved labeling in final printed form to the PMA Document Mail Center (HFZ-401) Center for Devices and Radiological Health, Food and Drug Administration (FDA), 9200 Corporate Boulevard, Rockville, Maryland 20850.

ADVERTISEMENT. No advertisement or other descriptive printed material issued by the applicant or private label distributor with respect to this device shall recommend or imply that the device may be used for any use that is not included in the FDA approved labeling for the device. If the FDA approval order has restricted the sale, distribution and use of the device to prescription use in accordance with 21 CFR 801.109 and specified that this restriction is being imposed in accordance with the provisions of section 520(e) of the act under the authority of section 515(d)(1)(B)(ii) of the act, all advertisements and other descriptive printed material issued by the applicant or distributor with respect to the device shall include a brief statement of the intended uses of the device and relevant warnings, precautions, side effects and contraindications.

PREMARKET APPROVAL APPLICATION (PMA) SUPPLEMENT. Before making any change affecting the safety or effectiveness of the device, submit a PMA supplement for review and approval by FDA unless the change is of a type for which a "Special PMA Supplement-Changes Being Effected" is permitted under 21 CFR 814.39(d) or an alternate submission is permitted in accordance with 21 CFR 814.39(e). A PMA supplement or alternate submission shall comply with applicable requirements under 21 CFR 814.39 of the final rule for Premarket Approval of Medical Devices.

All situations which require a PMA supplement cannot be briefly summarized, please consult the PMA regulation for further guidance. The guidance provided below is only for several key instances.

A PMA supplement must be submitted when unanticipated adverse effects, increases in the incidence of anticipated adverse effects, or device failures necessitate a labeling, manufacturing, or device modification.

A PMA supplement must be submitted if the device is to be modified and the modified device should be subjected to animal or laboratory or clinical testing designed to determine if the modified device remains safe and effective.

A "Special PMA Supplement - Changes Being Effectuated" is limited to the labeling, quality control and manufacturing process changes specified under 21 CFR 814.39(d)(2). It allows for the addition of, but not the replacement of previously approved, quality control specifications and test methods. These changes may be implemented before FDA approval upon acknowledgment by FDA that the submission is being processed as a "Special PMA Supplement by FDA Changes Being Effectuated." This acknowledgment is in addition to that issued by the PMA Document Mail Center for all PMA supplements submitted. This procedure is not applicable to changes in device design, composition, specifications, circuitry, software or energy source.

Alternate submissions permitted under 21 CFR 814.39(e) apply to changes that otherwise require approval of a PMA supplement before implementation of the change and include the use of a 30-day PMA supplement or annual postapproval report. FDA must have previously indicated in an advisory opinion to the affected industry or in correspondence with the applicant that the alternate submission is permitted for the change. Before such can occur, FDA and the PMA applicant(s) involved must agree upon any needed testing protocol, test results, reporting format, information to be reported, and the alternate submission to be used.

POSTAPPROVAL REPORTS. Continued approval of this PMA is contingent upon the submission of postapproval reports required under 21 CFR 814.84 at intervals of 1 year from the date of approval of the original PMA. Postapproval reports for supplements approved under the original PMA, if applicable, are to be included in the next and subsequent annual reports for the original PMA unless specified otherwise in the approval order for the PMA supplement. Two copies identified as "Annual Report" and bearing the applicable PMA reference number are to be submitted to the PMA Document Mail Center (HFZ-401), Center for Devices and Radiological Health, Food and Drug Administration, 9200 Corporate Boulevard, Rockville, Maryland 20850. The postapproval report shall indicate the beginning and ending date of the period covered by the report and shall include the following information required by 21 CFR 814.84:

- (1) Identification of changes described in 21 CFR 814.39(a) and changes required to be reported to FDA under 21 CFR 814.39(b).
- (2) Bibliography and summary of the following information not previously submitted as part of the PMA and that is

known to or reasonably should be known to the applicant:

- (a) unpublished reports of data from any clinical investigations or nonclinical laboratory studies involving the device or related devices ("related" devices include devices which are the same or substantially similar to the applicant's device); and
- (b) reports in the scientific literature concerning the device.

If, after reviewing the bibliography and summary, FDA concludes that agency review of one or more of the above reports is required, the applicant shall submit two copies of each identified report when so notified by FDA.

In addition to the above and in order to provide continued reasonable assurance of the safety and effectiveness of the device for its intended use, the annual postapproval reports shall include, separately for each model number (if applicable), the following information known by or reported to the applicant:

- (1) The number of pulse generators domestically implanted and the number of reported explants and deaths.
- (2) A breakdown of the reported deaths into pulse generator related and non-pulse generator related.
- (3) A breakdown of the reported explants into the numbers reported at end of battery life, having complications unresolvable by programming and for other reasons with safety and effectiveness issues which can be derived from the reports stated.
- (4) The number of pulse generators returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion and failed, with the failure mechanisms described.
- (5) A cumulative survival table for the pulse generators.
- (6) The number of programmers and modules shipped and the number of returns with a breakdown into the numbers currently in analysis, operating properly and failed, with the failure mechanisms described.

ADVERSE REACTION AND DEVICE DEFECT REPORTING. As provided by 21 CFR 814.82(a)(9), FDA has determined that in order to provide continued reasonable assurance of the safety and effectiveness of the device, the applicant shall submit 3 copies of a written report identified, as applicable, as an "Adverse Reaction Report" or "Device Defect Report" to the PMA Document Mail Center (HFZ-401), Center for Devices and Radiological Health, Food and Drug Administration, 9200 Corporate Boulevard, Rockville, Maryland 20850 within 10 days after the applicant receives or has knowledge of information concerning:

- (1) A mix-up of the device or its labeling with another article.
- (2) Any adverse reaction, side effect, injury, toxicity, or sensitivity reaction that is attributable to the device and
 - (a) has not been addressed by the device's labeling or
 - (b) has been addressed by the device's labeling, but is occurring with unexpected severity or frequency.
- (3) Any significant chemical, physical or other change or deterioration in the device or any failure of the device to meet the specifications established in the approved PMA that could not cause or contribute to death or serious injury but are not correctable by adjustments or other maintenance procedures described in the approved labeling. The report shall include a discussion of the applicant's assessment of the change, deterioration or failure and any proposed or implemented corrective action by the applicant. When such events are correctable by adjustments or other maintenance procedures described in the approved labeling, all such events known to the applicant shall be included in the Annual Report described under "Postapproval Reports" above unless specified otherwise in the conditions of approval to this PMA. This postapproval report shall appropriately categorize these events and include the number of reported and otherwise known instances of each category during the reporting period. Additional information regarding the events discussed above shall be submitted by the applicant when determined by FDA to be necessary to provide continued reasonable assurance of the safety and effectiveness of the device for its intended use.

REPORTING UNDER THE MEDICAL DEVICE REPORTING (MDR) REGULATION.

The Medical Device Reporting (MDR) Regulation became effective on December 13, 1984. This regulation was replaced by the reporting requirements of the Safe Medical Devices Act of 1990 which became effective July 31, 1996, and requires that all manufacturers and

importers of medical devices, including in vitro diagnostic devices, report to the FDA whenever they receive or otherwise become aware of information, from any source, that reasonably suggests that a device marketed by the manufacturer or importer:

- (1) May have caused or contributed to a death or serious injury; or
- (2) Has malfunctioned and such device or similar device marketed by the manufacturer or importer would be likely to cause or contribute to a death or serious injury if the malfunction were to recur.

The same events subject to reporting under the MDR Regulation may also be subject to the above "Adverse Reaction and Device Defect Reporting" requirements in the "Conditions of Approval" for this PMA. FDA has determined that such duplicative reporting is unnecessary. Whenever an event involving a device is subject to reporting under both the MDR Regulation and the "Conditions of Approval" for a PMA, the manufacturer shall submit the appropriate reports required by the MDR Regulation within the time frames as identified in 21 CFR 803.10(c) using FDA Form 3500A, i.e., 30 days after becoming aware of a reportable death, serious injury, or malfunction as described in 21 CFR 803.50 and 21 CFR 803.52 and 5 days after becoming aware that a reportable MDR event requires remedial action to prevent an unreasonable risk of substantial harm to the public health. The manufacturer is responsible for submitting a baseline report on FDA Form 3417 for a device when the device model is first reported under 21 CFR 803.50. This baseline report is to include the PMA reference number. Any written report and its envelope is to be specifically identified, e.g., "Manufacturer Report," "5-Day Report," "Baseline Report," etc. Any written report is to be submitted to:

Food and Drug Administration
Center for Devices and Radiological Health
Medical Device Reporting
PO Box 3002
Rockville, Maryland 20847-3002

Copies of the MDR Regulation (FOD # 336&1336) and FDA publications entitled "An Overview of the Medical Device Reporting Regulation" (FOD # 509) and "Medical Device Reporting for Manufacturer" (FOD # 987) are available on the CDRH WWW Home Page. They are also available through CDRH's Fact-On-Demand (F-O-D) at 800-899-0381. Written requests for information can be made by sending a facsimile to CDRH's Division of Small Manufacturers Assistance (DSMA) at 301-443-8818.

APPENDIX B: SYSTEM EXPERIENCE REPORT

**VENTAK AV/VR PRIZM
PULSE GENERATORS
EXTERNAL DEVICES**

PMA NO. P960040

Cardiac Pacemakers, Inc.
4100 Hamline Avenue North
St. Paul, MN 55112

Submission Date: August 13, 2001

Prepared and Submitted by:

Lynn C. Zurn

Name

Signature

Reliability Assurance Event Analyst

Title

Date

8-13-01

Armando Beiro

Name

Signature

Reliability Assurance Engineer

Title

Date

8-14-01

Dan Tich

Name

Signature

Reliability Assurance Manager

Title

Date

8/13/01

TABLE OF CONTENTS

I. REPORT OVERVIEW	3
SCOPE OF REPORT	3
II. DEVICE EXPERIENCE	4
VENTAK AV, AICD PULSE GENERATOR, MODEL 1810	5
VENTAK AV AICD PULSE GENERATOR, MODEL 1720	9
VENTAK AV II DDD AICD PULSE GENERATOR, MODEL 1820	12
VENTAK AV II DR AICD PULSE GENERATOR, MODEL 1821	15
VENTAK AV II DR AICD PULSE GENERATOR, MODEL 1826	19
VENTAK AV III DR AICD PULSE GENERATOR, MODEL 1831	22
VENTAK AV III DR AICD, PULSE GENERATOR, MODEL 1836	28
VENTAK VR 3.2 AICD PULSE GENERATOR, MODEL 1774	32
VENTAK VR 6.1 AICD PULSE GENERATOR, MODEL 1775	36
VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1850	39
VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1851	43
VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1856	48
VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1853	51
VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1858	55
VENTAK PRIZM VR/HE PULSE GENERATOR, MODEL 1852	58
VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1857	61
VENTAK PRIZM 2 DR/VR PULSE GENERATOR, MODEL 1861	64
II. EXTERNALS EXPERIENCE	67
SOFTWARE DISK FOR THE VENTAK AV & AV II DDD PULSE GENERATORS, MODEL 2833	68
SOFTWARE DISK FOR THE VENTAK AV II & III DR PULSE GENERATORS, MODEL 2843	68
SOFTWARE DISK FOR THE VENTAK VR PULS E GENERATORS, MODEL 2841	68
SOFTWARE DISK FOR THE VENTAK PRIZM DR/VR PULSE GENERATORS, MODEL 2844	69

I. REPORT OVERVIEW

SCOPE OF REPORT

This Ventak AV/VR PRIZM Experience System Report presents performance information for the Ventak AV/VR PRIZM Pulse Generator Models and External Devices approved under PMA P960040. Device data is segregated by model number to facilitate review.

This report includes device performance information for the period between June 1, 2000 through May 31, 2001. The cumulative survival analysis is an ongoing accumulation of data on device life since the first documented implant.

<u>Model</u>	<u>Device Name</u>
1810	Ventak AV, AICD Pulse Generator
1815	Ventak AV, AICD Pulse Generator
1820	Ventak AV II DDD, AICD Pulse Generator
1825*	Ventak AV II DDD, AICD Pulse Generator
1821	Ventak AV II DR, AICD Pulse Generator
1826	Ventak AV II DR, AICD Pulse Generator
1831	Ventak AV III DR, AICD Pulse Generator
1836	Ventak AV III DR, AICD Pulse Generator
1774	Ventak VR, AICD Pulse Generator
1775	Ventak VR, AICD Pulse Generator
1850	Ventak PRIZM DR/VR, Pulse Generator
1851	Ventak PRIZM DR/VR, Pulse Generator
1852	Ventak PRIZM DR/VR, Pulse Generator
1855*	Ventak PRIZM DR/VR, Pulse Generator
1856	Ventak PRIZM DR/VR, Pulse Generator
1853	Ventak PRIZM DR/HE, Pulse Generator
1858	Ventak PRIZM DR/HE, Pulse Generator
1852	Ventak PRIZM VR/HE, Pulse Generator
1857	Ventak PRIZM VR/HE, Pulse Generator
1861	Ventak PRIZM 2 DR/VR, Pulse Generator
1862*	Ventak PRIZM 2 DR/VR, Pulse Generator
2833	Software Disk for Ventak AV & AV II DDD Pulse Generators
2843	Software Disk for Ventak AV II & III DR Pulse Generators
2841	Software Disk for Ventak VR Pulse Generators
2844	Software Disk for Ventak PRIZM DR/VR Pulse Generators

**(These models had no activity during this report period.)*

II. DEVICE EXPERIENCE

The following information is submitted in response to requirements as specified in the Conditions of Approval. The numbers reported here are for devices distributed in the United States. Data used to derive these reports were obtained from the company's databases. All events tabulated in this report were reviewed for Medical Device Reporting (MDR) filing at the time GUIDANT received information that a reportable event may have occurred. Where appropriate, MDRs were filed.

Data are reported separately for each model as specified in the Conditions of Approval.

The cumulative survival tables in this report were prepared using standard life table techniques (see Cutler and Ederer, 1958, J Chron. Dis., 8:6). The follow-up experience has been divided into three-month intervals for reporting purposes. For each interval, a calculation has been made of the number of units at risk during the interval. The number of units at risk was found by starting with the number of devices that entered the interval and correcting for the units that failed, or were taken out-of-service for a reason not related to the device.

The number of devices that failed during an interval, divided by the number at risk, result in an estimate of the probability of failure during the interval.

Please note that the Cumulative Survival Percentage provided in this report is affected only by those units that were returned to GUIDANT, analyzed, determined to be defective and classified as such by Reliability Assurance. Units with induced damage (such as "terminal pin bent"), that are considered out of specification even though they are operating within specification, do not factor into the reported Cumulative Survival Percent.

For each interval, the probability of failure during the interval is subtracted from one to give the probability of survival in the interval. Survival probability through a series of intervals is obtained as the product of individual survival probabilities.

Please note, that the total number of units reported on the cumulative survival report table may be higher than the total number of units reported under FDA Reporting Requirement #1. The report period for FDA Reporting Requirement #1 represents only a portion of the total number of units implanted for that model, while cumulative survival represents all devices implanted since release for clinical use.

VENTAK AV, AICD PULSE GENERATOR, MODEL 1810

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	91
(c)	Total units in patients who reportedly died in report period	28

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device related	28
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	91
Elective Replacement	24
Erosion	1
Heart Transplant	2
Infection	2
Normal ERI	19
Other observation/Complication	1
Product Performance Issue**	9
Communication/Telemetry Issue	3
Diagnostic/Data Issue	1
Sensing/Detection Issue	5
System to Patient Interface Issue	3
Therapy Delivery/Effectiveness Issue	4
Unknown or Insufficient information provided	31
Upgrade	2

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	43
Out of Specification	6
Arcing damage, Module output high power	1
Induced	2
Overstress, high energy	1
Solder joint, cracked	1
Undetermined/inconclusive	1
Unit Meet Specifications	37

(b) Analysis results of preimplant units returned for cause	0
---	---

5. FDA Reporting Requirement: A cumulative survival table for the devices.

The Cumulative Survival Analysis for VENTAK AV AICD Pulse Generator, Model 1810

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		99.82%	0.07%
3 - 6		99.77%	0.08%
6 - 9		99.71%	0.09%
9 - 12		99.60%	0.10%
12 - 15		99.57%	0.11%
15 - 18		99.51%	0.12%
18 - 21		99.48%	0.12%
21 - 24		99.39%	0.13%
24 - 27		99.36%	0.13%
27 - 30		99.21%	0.15%
30 - 33		99.14%	0.16%
33 - 36		99.08%	0.16%
36 - 39		99.01%	0.17%
39 - 42		99.01%	0.17%
42 - 45		98.95%	0.18%
45 - 48		98.95%	0.18%

VENTAK AV AICD PULSE GENERATOR, MODEL 1720

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	6
(c)	Total units in patients who reportedly died in report period	1

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the device may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device related	1
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	6
Heart Transplant	1
Normal ERI	2
Unknown of insufficient information provided	3

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as leads not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	3
Unit Meets Specification	3

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for VENTAK AV AICD Pulse Generator, Model 1815

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	[REDACTED]	100.00%	0.00%
3 - 6	[REDACTED]	100.00%	0.00%
6 - 9	[REDACTED]	100.00%	0.00%
9 - 12	[REDACTED]	100.00%	0.00%
12 - 15	[REDACTED]	100.00%	0.00%
15 - 18	[REDACTED]	99.36%	0.63%
18 - 21	[REDACTED]	99.36%	0.63%
21 - 24	[REDACTED]	98.70%	0.91%
24 - 27	[REDACTED]	98.70%	0.91%
27 - 30	[REDACTED]	98.70%	0.91%
30 - 33	[REDACTED]	98.70%	0.91%
33 - 36	[REDACTED]	98.70%	0.91%
36 - 39	[REDACTED]	98.70%	0.91%
39 - 42	[REDACTED]	98.70%	0.91%
42 - 45	[REDACTED]	98.70%	0.91%

VENTAK AV II DDD AICD PULSE GENERATOR, MODEL 1820

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	4
(c)	Total units in patients who reportedly died in report period	5

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device related	5
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	4
Elective Replacement	1
Normal ERI	1
Unknown of insufficient information provided	2

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as leads not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	4
Out of Specification	1
Cracked, structure or outer body, Analog hybrid Assy 1	
Unit Meet Specifications	3
(b) Analysis results of preimplant units returned for cause	0

5. FDA Reporting Requirement: A cumulative survival table for the devices.

The Cumulative Survival Analysis for VENTAK AV II DDD AICD Pulse Generator, Model 1820

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	████████	100.00%	0.00%
3 - 6	████████	99.82%	0.18%
6 - 9	████████	99.82%	0.18%
9 - 12	████████	99.82%	0.18%
12 - 15	████████	99.82%	0.18%
15 - 18	████████	99.82%	0.18%
18 - 21	████████	99.82%	0.18%
21 - 24	████████	99.82%	0.18%
24 - 27	████████	99.62%	0.27%
27 - 30	████████	99.62%	0.27%
30 - 33	████████	99.62%	0.27%
33 - 36	████████	99.62%	0.27%
36 - 39	████████	99.62%	0.27%
39 - 42	████████	99.62%	0.27%

VENTAK AV II DR AICD PULSE GENERATOR, MODEL 1821

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	83
(c)	Total units in patients who reportedly died in report period	36

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the device may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	36
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	83
Dissatisfied with product	1
Coincident Removal	3
Elective Replacement	27
Heart Transplant	3
Infection	2
Normal ERI	7
Product Performance Issue**	10
Communication/Telemetry Issue	2
Diagnostic/ Data Issue	4
Mechanical Connection Issue	1
Premature Battery Depletion	2
Sensing/ Detection Issue	2
System to Patient Interface Issue	2
Therapy Delivery/ Effectiveness Issue	2
Unknown or Insufficient Information provided	28
Upgrade	2

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as leads not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause		52
Out of Specification		16
Transfer Error Code	1	
Cracked substrate hybrid assy Analog AVII	1	
Depletion, premature, undetermined	1	
Induced	8	
Logic Error	1	
Overstress, high energy	1	
Shorted, module output high power	1	
Solder Joint, poor quality	1	
Undetermined/inconclusive	1	
Unit Meet Specifications		36

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for VENTAK AV II DR AICD Pulse Generator, Model 1821

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		99.96%	0.03%
3 - 6		99.84%	0.06%
6 - 9		99.60%	0.09%
9 - 12		99.58%	0.09%
12 - 15		99.54%	0.10%
15 - 18		99.54%	0.10%
18 - 21		99.52%	0.10%
21 - 24		99.50%	0.10%
24 - 27		99.50%	0.10%
27 - 30		99.45%	0.11%
30 - 33		99.40%	0.11%
33 - 36		99.29%	0.13%
36 - 39		99.29%	0.13%

VENTAK AV II DR AICD PULSE GENERATOR, MODEL 1826

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	7
(c)	Total units in patients who reportedly died in report period	1

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the device may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	1
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	7
Normal ERI	2
Unknown of insufficient information provided	5

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	0
--	---

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for VENTAK AV II DR AICD Pulse Generator, Model 1826

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	[REDACTED]	100.00%	0.00%
3 - 6	[REDACTED]	100.00%	0.00%
6 - 9	[REDACTED]	100.00%	0.00%
9 - 12	[REDACTED]	100.00%	0.00%
12 - 15	[REDACTED]	100.00%	0.00%
15 - 18	[REDACTED]	100.00%	0.00%
18 - 21	[REDACTED]	100.00%	0.00%
21 - 24	[REDACTED]	100.00%	0.00%
24 - 27	[REDACTED]	100.00%	0.00%
27 - 30	[REDACTED]	100.00%	0.00%
30 - 33	[REDACTED]	100.00%	0.00%
33 - 36	[REDACTED]	100.00%	0.00%



VENTAK AV III DR AICD PULSE GENERATOR, MODEL 1831

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	██████
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	188
(c)	Total units in patients who reportedly died in report period	84

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the device may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	84
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

See table of following page for details

Out-of-Service Units (as indicated by the user)	188
Coincident Removal	2
Elective Replacement	21
Infection	12
Normal ERI	2
Other Observation/ complication	1
Heart Transplant	14
Product Performance Issue**	78
Communication/Telemetry Issue	40
Contamination Issue	1
Diagnostic/Data Issue	57
Mechanical Connection Issue	1
Premature Battery Depletion	14
Sensing/ Detection Issue	4
System to Patient Interface Issue	8
Therapy Delivery/Effectiveness Issue	21
Product Performance Issue, Never Implanted	1
Unknown or Insufficient information provided	56
Upgrade	1

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	144
Analysis Pending	18

Out of Specification	67
Degraded, parameter shift out of spec	3
Dendritic growth	1
Dendritic growth, EPROM	1
Dendritic growth, eprom encoded	2
Dendritic growth, hybrid assembly	12
Depletion, premature, undetermined	1
Excessive current, undetermined	1
Fracture, non specific, inductor power	1
Fracture, non specific, LDFrame	1
Induced	9
Leaky, Cap Tant Module Dual	2
Leaky, Cap, Tantalum	2
Leaky Hybrid assembly	1
Leaky, Trans Die Mosfet NCHAN Sense	1
Logic error	1
Logic error, Firmware/Software	3
Memory/address value(s) corrupted	1
Severed, LD Frame/ASIC Telem	1
Shorted, dendritic growth, Eprom Encoded	1
Shorted, dendritic growth, hybrid assembly	1
Shorted, foreign material, XFMR Toroid	1
Solder joint, electrically intermittent	2
Undetermined/inconclusive	21
Undetermined/inconclusive Hybrid assembly	1
Unit Meet Specifications	59

(b) Analysis results of preimplant units returned for cause	3
Out of Specification	1
Induced	1
Unit Meet Specifications	2

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak AV III DR AICD Pulse Generator, Model 1831

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		99.90%	0.03%
3 - 6		99.80%	0.05%
6 - 9		99.63%	0.06%
9 - 12		99.38%	0.08%
12 - 15		99.16%	0.09%
15 - 18		99.03%	0.10%
18 - 21		98.86%	0.11%
21 - 24		98.77%	0.12%
24 - 27		98.58%	0.14%
27 - 30		98.42%	0.16%

VENTAK AV III DR AICD, PULSE GENERATOR, MODEL 1836

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	██████
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	12
(c)	Total units in patients who reportedly died in report period	4

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its devices may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the device may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	4
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	12
Elective Replacement	4
Infection	1
Product Performance Issue**	1
Communication/ Telemetry Issue	1
System to Patient Interface Issue	1
Unknown or insufficient information provided	5
Upgrade	1

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) not every device taken out-of-service is returned for analysis, and 2) pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	5
Analysis Pending	1
Unit Meet Specifications	4

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak AV III DR AICD Pulse Generator, Model 1836

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		100.00%	0.00%
3 - 6		100.00%	0.00%
6 - 9		99.22%	0.78%
9 - 12		99.22%	0.78%
12 - 15		99.22%	0.78%
15 - 18		99.22%	0.78%
18 - 21		99.22%	0.78%
21 - 24		99.22%	0.78%
24 - 27		99.22%	0.78%

VENTAK VR 3.2 AICD PULSE GENERATOR, MODEL 1774

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	16
(c)	Total units in patients who reportedly died in report period	8

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	8
• Device Related	0

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4

(a) Analysis results of explanted units returned for cause	12
Out of Specification	3
Arcing damage, XFMR Toroid	1
Induced	1
Undetermined/inconclusive	1
Unit Meet Specifications	9

(b) Analysis results of preimplant units returned for cause	1
Unit Meet Specifications	1

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for VENTAK VR 3.2 AICD Pulse Generator, Model 1774

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	[REDACTED]	99.89%	0.11%
3 - 6	[REDACTED]	99.77%	0.16%
6 - 9	[REDACTED]	99.65%	0.20%
9 - 12	[REDACTED]	99.65%	0.20%
12 - 15	[REDACTED]	99.65%	0.20%
15 - 18	[REDACTED]	99.65%	0.20%
18 - 21	[REDACTED]	99.65%	0.20%
21 - 24	[REDACTED]	99.65%	0.20%

OK

VENTAK VR 6.1 AICD PULSE GENERATOR, MODEL 1775

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	██████
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	7
(c)	Total units in patients who reportedly died in report period	1

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	1
• Device Related	0

68

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	7
Elective Replacement	3
Unknown or insufficient information provided	4

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	6
Unit Meet Specifications	6

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for VENTAK VR 6.1 AICD Pulse Generator, Model 1775

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		100.00%	0.00%
3 - 6		100.00%	0.00%
6 - 9		100.00%	0.00%
9 - 12		100.00%	0.00%
12 - 15		100.00%	0.00%
15 - 18		100.00%	0.00%
18 - 21		100.00%	0.00%
21 - 24		100.00%	0.00%

VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1850

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	135
(c)	Total units in patients who reportedly died in report period	126

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	126
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

See Table on Following page for details

Out-of-Service Units (as indicated by the user)	135
Elective Replacement	26
Electively not Used	1
Erosion	2
Heart Transplant	9
Infection	21
Normal ERI	2
Other Observation/Complication	3
Product Performance Issue**	22
Communication/Telemetry Issue	3
Diagnostic/Data Issue	19
Mechanical Connection Issue	2
Premature Battery Depletion	4
Sensing/ Detection Issue	1
System to Patient Interface Issue	9
Therapy Delivery/ Effectiveness Issue	3
Recall/Advisory	1
Unknown or insufficient information provided	44
Upgrade	4

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

72

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	65
Analysis Pending	1
Out of Specification	29
Degraded, parameter shift, Hybrid	1
Degraded, parameter shift, Transistor	1
Induced	5
Leaky, Cap Chip Tant, Module 5	1
Logic Error	1
Memory/address value(s) corrupted	18
Trace or via, open, Hybrid MDL	1
Undetermined/inconclusive	1
Unit Meet Specifications	35

ing

(b) Analysis results of preimplant units returned for cause	91
Analysis Pending	2
Out of Specification	22
Degraded, parameter shift, Trans-Die Mosfet	1
Depletion, premature, undetermined	1
Induced	11
Memory/address value(s) corrupted	5
Trace or via, open, hybrid MDL	2
Undetermined/inconclusive	2
Unit Meet Specifications	67

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM DR/VR Pulse Generator, Model 1850

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	[REDACTED]	99.83%	0.05%
3 - 6	[REDACTED]	99.59%	0.09%
6 - 9	[REDACTED]	99.53%	0.10%
9 - 12	[REDACTED]	99.53%	0.10%
12 - 15	[REDACTED]	99.53%	0.10%

VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1851

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	180
(c)	Total units in patients who reportedly died in report period	294

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	293
• Device Related (1851/309036) see below FDA # 2124215000200002428	1

Guidant received information that the patient with this Implantable Cardioverter Defibrillator (ICD) expired 2 days following the implant of this device. It was noted that the device "went off" continuously for approximately 10 minutes; however, it was not noted as to whether or not this occurred prior to or after the patient's death. Guidant made several attempts to obtain additional information regarding the exact cause of the patient's death and the performance of the device; however, additional information was unable to be obtained. The device was not returned for analysis as it was buried with the patient. Guidant will submit additional information if it becomes available. An MDR death report was submitted November 1, 2000.

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	180
Elective Replacement	17
Heart Transplant	13
Infection	47
Normal ERI	3
Other Observation/Complication	4
Product Performance Issue**	37
Communication/Telemetry Issue	1
Diagnostic/ Data Issue	40
Mechanical Connection Issue	17
Physical Damage	1
Premature Battery Depletion	2
Sensing/Detection Issue	6
System to Patient Interface Issue	15
Therapy Delivery/Effectiveness Issue	4
Recall/Advisory	2
Unknown or insufficient information provided	55
Upgrade	2

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause		109
Analysis Pending		1
Out of Specification		55
Degraded, parameter shift out of spec	1	
Foreign Material present, header ICD	1	
Induced	15	
Leaky, Cap Tant Module dual	1	
Logic error	1	
Memory/address value(s) corrupted	24	
Open, electrically, hybrid MDL	2	
Shorted, contact with case	2	
Solder joint, cracked	1	
Solder joint, electrically intermittent	1	
Trace or via open, Hybrid MDL	6	
Unit Meet Specifications		53

(b) Analysis results of preimplant units returned for cause		107
Analysis Pending		3
Out of Specification		39
Foreign material present, header ICD	1	
Induced	23	
Memory/address value(s) corrupted	9	
Open, electrically	1	
Stuck, set screw socket head half dog	1	
Trace or via, open, hybrid MDL	2	
Undetermined/inconclusive	2	
Unit Meet Specifications		65

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM HE VR Pulse Generator, Model 1851

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		99.87%	0.04%
3 - 6		99.72%	0.05%
6 - 9		99.66%	0.06%
9 - 12		99.61%	0.07%
12 - 15		99.61%	0.07%

VENTAK PRIZM DR/VR PULSE GENERATOR, MODEL 1856

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	0
(c)	Total units in patients who reportedly died in report period	0

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	0
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	0
---	---

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	0
--	---

(b) Analysis results of preimplant units returned for cause	0
---	---

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

**The Cumulative Survival Analysis for Ventak PRIZM Pulse Generator,
Model 1856**

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3		100.00%	0.00%
3 - 6		100.00%	0.00%
6 - 9		100.00%	0.00%
9 - 12		100.00%	0.00%

VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1853

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	18
(c)	Total units in patients who reportedly died in report period	34

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	34
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	18
Heart Transplant	6
Infection	6
Product Performance Issue**	2
Diagnostic/ Data Issue	3
Mechanical Connection Issue	3
System to Patient Interface Issue	2
Recall/Advisory	1
Unknown or Insufficient information provided	3

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	12
Out of Specification	5
Induced 4	
Trace or via, open hybrid MDL 1	
Unit Meet Specifications	7

(b) Analysis results of preimplant units returned for cause	7
Out of Specification	3
Induced 2	
Memory/address value(s) corrupted 1	
Unit Meet Specifications	4

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM DR/HE Pulse Generator, Model 1853

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0-3	██████████	99.89%	0.11%
3-6	██████████	89.89%	0.11%
6-9	██████████	99.89%	0.11%

VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1858

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

* GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a) Total units implanted in report period	
(b) Total units reported explanted (out-of-service) in report period (excluding reported deaths)	1
(c) Total units in patients who reportedly died in report period	1

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	1
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	1
Elective Replacement	1

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	0
(b) Analysis results of preimplant units returned for cause	4
Out of Specification	1
Foreign material present	1
Unit Meet Specification	3

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM DR/HE Pulse Generator, Model 1858

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0-3	██████	100.00%	0.00%
3-6	██████	100.00%	0.00%
6-9	██████	100.00%	0.00%

VENTAK PRIZM VR/HE PULSE GENERATOR, MODEL 1852

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	5
(c)	Total units in patients who reportedly died in report period	2

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	2
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	5
Heart Transplant	1
Product Performance Issue**	2
Diagnostic/Data Issue	3
Unknown or Insufficient information provided	2

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	1
Out of Specification	1
Memory/address value(s) corrupted	1

(b) Analysis results of preimplant units returned for cause	0
---	---

at

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM VR/HE Pulse Generator, Model 1852

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0-3	[REDACTED]	100.00%	0.00%
3-6	[REDACTED]	99.21%	0.79%

VENTAK PRIZM DR/HE PULSE GENERATOR, MODEL 1857

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	1
(c)	Total units in patients who reportedly died in report period	3

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	3
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	1
Elective Replacement	1

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	1
Unit Meet Specification	1
(b) Analysis results of preimplant units returned for cause	1
Out of Specification	1
Memory/address value(s) corrupted	1

99

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM VR/HE Pulse Generator, Model 1857

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0-3		100.00%	0.00%
3-6		100.00%	0.00%

VENTAK PRIZM 2 DR/VR PULSE GENERATOR, MODEL 1861

1. *FDA Reporting Requirement: The number of devices domestically implanted and the number of reported explants* and deaths.*

*GUIDANT uses the term "out-of-service" rather than "explanted" because GUIDANT does not always know if the device was removed from the body.

(a)	Total units implanted in report period	
(b)	Total units reported explanted (out-of-service) in report period (excluding reported deaths)	50
(c)	Total units in patients who reportedly died in report period	79

2. *FDA Reporting Requirement: A breakdown of the reported deaths into device related and non-device related.*

GUIDANT used the company's Medical Device Reporting procedure to determine death events that met the definition of "Device Related."

A death is reported as "Device Related" when GUIDANT received or otherwise became aware of information that reasonably suggested that one of its leads may have caused or contributed to the death. In addition, a death is reported as "Device Related" if the death was believed to be the result of user error when the information reasonably suggested that the lead may have caused or contributed to the death.

Submission of a report as "Device Related" does not constitute an admission that the manufacturer, product, medical personnel, or user facility caused or contributed to the death.

• Non-Device Related	79
• Device Related	0

3. *FDA Reporting Requirement: A breakdown of the reported explants into the numbers reported as at end of battery life, having complications unresolvable by programming, and for other reasons with safety and effectiveness issues which can be derived from the reports stated.*

Please note that GUIDANT designates units as "out-of-service" rather than "explanted." When an out-of-service device is explanted and returned to GUIDANT, it is tabulated in the analysis section (FDA Reporting Requirement #4).

Out-of-Service Units (as indicated by the user)	50
Elective Replacement	4
Electively not Used	3
Heart Transplant	1
Infection	12
Other observation/complication	4
Product Performance Issue**	11
Communication/Telemetry Issue	1
Diagnostic/ Data Issue	6
Mechanical Connection Issue	5
Premature Battery Depletion	1
Sensing/ Detection Issue	3
System to patient Interface Issue	11
Therapy Delivery/ Effectiveness Issue	2
Unknown or Insufficient information provided	15

***Each Out-of-Service unit with a product performance issue may have more than one observation. Therefore, the number of units with product performance issues may be less than the number of observations.*

4. *FDA Reporting Requirement: The number of devices returned to the applicant for cause from domestic sources with a breakdown into the numbers currently in analysis, operating properly, at normal battery depletion, and failed, with the failure mechanisms described.*

Please note that in this section GUIDANT is reporting devices that have been explanted and returned for analysis, as well as devices not implanted before return. Devices in the later category are considered pre-implant units and are reported separately.

The number of devices discussed in this section will differ from those discussed in FDA Reporting Requirement #3 for two reasons: 1) Not every device taken out-of-service is returned for analysis, and 2) Pre-implant devices are included in FDA Reporting Requirement #4.

(a) Analysis results of explanted units returned for cause	22
Analysis Pending	3
Out of Specification	6
Cracked, RC Network	1
Degraded, parameter shift out of spec	1
Induced	3
Leaky, Cap Tantalum	1
Unit Meet Specifications	13

(b) Analysis results of preimplant units returned for cause	58
Analysis Pending	6
Out of Specification	8
Induced	8
Unit Meet Specification	44

5. *FDA Reporting Requirement: A cumulative survival table for the devices.*

The Cumulative Survival Analysis for Ventak PRIZM 2 DR/VR Pulse Generator, Model 1861

AGE in MONTHS	NUMBER of UNITS	SURVIVAL RATE	STANDARD ERROR
0 - 3	[REDACTED]	99.98%	0.02%
3 - 6	[REDACTED]	99.88%	0.07%

EXTERNAL PRODUCTS EXPERIENCE

I.FDA Reporting Requirement: The number of programmers and modules shipped and the number of returns with a breakdown in to the numbers currently in analysis, operating properly and failed, with the failure mechanisms described.

SOFTWARE DISK FOR THE VENTAK AV & AV II DDD PULSE GENERATORS, MODEL 2833

(a)	Total units shipped in report period	
(b)	Analysis results of explanted units returned for cause	101
	Scrap	80
	Unit Meet Specifications	21

SOFTWARE DISK FOR THE VENTAK AV II & III DR PULSE GENERATORS, MODEL 2843

(a)	Total units shipped in report period	
(b)	Analysis results of total units returned in report period	88
	Scrap	88

SOFTWARE DISK FOR THE VENTAK VR PULS E GENERATORS, MODEL 2841

(a)	Total units shipped in report period	
(b)	Analysis results of total units returned in report period	107
	Scrap	81
	Unit Meets Specifications	26

**SOFTWARE DISK FOR THE VENTAK PRIZM DR/VR PULSE GENERATORS,
MODEL 2844**

(a)	Total units shipped in report period	
(b)	Analysis results of total units returned in report period	236
	Scrap	194
	Analysis Pending	2
	Out of Specification	1
	Degraded, parameter out of spec	1
	Unit Meet Specification	39

701

APPENDIX C: CAPACITOR LIFE TEST REPORT for VENTAK VR

102

GUIDANT

[REDACTED], Revision C

**Component Evaluation Test Report
High Energy Capacitor
Storage Test- Interim Test Results
CPI Part Number [REDACTED]**

Manufacturer: Guidant
Project: Ventak VR
Responsible Department: Supplier Development

Authored By:

Maureen Esteban

(name)

(signature)

Component Engineer

(title)

7/10/01

(date)

Approved By:

Jayne Fangel

(name)

(signature)

Supplier Development Mgr.

(title)

7/19/01

(date)

103

21 Whole Page Redactions