


CLIENT LOGIN

Site Test Catalog

SEARCH



- Online Test Catalog Home
- New Tests
- Test Updates
- Special Request Testing
- Innovation
- 40th Anniversary Home
- Suggest a Test
- Tox Wiki
- Blog
- Site Map

Ephedrines Panel, Blood Test (4023B)

Analysis Code 4023B

Test Name Ephedrines Panel, Blood

Test Includes Ephedrine [LC-MS/MS], Methylephedrine [LC-MS/MS], Norpseudoephedrine [LC-MS/MS], Phenylpropanolamine [LC-MS/MS], Pseudoephedrine [LC-MS/MS]

Compound Synonym(s) Cathine Ephedra Ma Huang Norephedrine PPA PSE Sudafed®

Purpose Identification and Quantitation by LC-MS/MS

Category Antihistamine, Decongestant, Stimulant, Anorexogenic, Bronchodilator, Stimulant, Decongestant, Stimulant

Method(s) High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Specimen Requirements 1 mL Blood

Transport Temperature Refrigerated

Specimen Container NMS Labs has no experimental or literature-based data regarding the choice of specific specimen collection containers for this test.

Special Handling None

Light Protection Required Not Required

Stability Room Temperature: 16 day(s)
Refrigerated: 16 day(s) Frozen (-20 °C): 16 day(s)

Rejection Criteria None

Day(s) Test Set-up [LC-MS/MS] Monday Wednesday Friday
2nd Shift 3 days

Suggested CPT Code [LC-MS/MS] 83789

View and compare Ephedrines Panel, Blood Test (4023B)


ShareThis

A	B	C	D	E	F	G	H	I	J	K	L	M
N	O	P	Q	R	S	T	U	V	W	X	Y	Z
0	1	2	3	4	5	6	7	8	9			
<input type="button" value="SEARCH"/>												

SIGN ME UP




Sign me up for e-mail updates



[CLIENT LOGIN](#)

[Site](#) [Test Catalog](#)

[SEARCH](#)



- [Online Test Catalog Home](#)
- [New Tests](#)
- [Test Updates](#)
- [Special Request Testing](#)
- [Innovation](#)
- [40th Anniversary Home](#)
- [Suggest a Test](#)
- [Tox Wiki](#)
- [Blog](#)
- [Site Map](#)

Amphetamines Panel, Blood (Forensic) Test (8215B)

Analysis Code 8215B

Test Name Amphetamines Panel, Blood (Forensic)

Test Includes Amphetamine [LC-MS/MS], Amphetamines [ELISA], Ephedrine [LC-MS/MS], MDA [LC-MS/MS], MDEA [LC-MS/MS], MDMA [LC-MS/MS], Methamphetamine [LC-MS/MS], Methylephedrine [LC-MS/MS], Norpseudoephedrine [LC-MS/MS], Phenmetrazine [LC-MS/MS], Phenmetrazine [LC-MS/MS], Phentermine [LC-MS/MS], Phenylpropanolamine [LC-MS/MS], Pseudoephedrine [LC-MS/MS], Selegiline [LC-MS/MS]

Compound Synonym(s) Adam Adipex-P® Bontril® Cathine Crank Crystal Deprenyl Desoxyn® Dexedrine Ecstasy Eldepryl® Ephedra Eve Ionamin® Ma Huang Methedrine Methylenedioxyamphetamine Methylenedioxyethylamphetamine Methylenedioxymethamphetamine Norephedrine PPA PSE Plegine® Prelu-2® Preludin® Pro-Fast® Speed Sudafed® Teramine® XTC

Purpose Forensic Analysis

Category Antihistamine, Decongestant, Stimulant, Stimulant, Anorexogenic, Appetite Suppressant, Bronchodilator, Stimulant, Decongestant, Stimulant

Method(s) Enzyme-Linked Immunosorbent Assay (ELISA) High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Specimen Requirements 2 mL Blood

Transport Temperature Refrigerated

Specimen Container NMS Labs has no experimental or literature-based data regarding the choice of specific specimen collection containers for this test.

Special Handling None

Light Protection Required Not Required

Stability Room Temperature: 2 day(s)
Refrigerated: 14 day(s) Frozen (-20 °C): 16 day(s)

Rejection Criteria Received Room Temperature.

Day(s) Test Set-up [ELISA] Monday-Saturday 2nd Shift 1 day
[LC-MS/MS] Monday-Friday 2nd Shift 3 days

Suggested CPT Code [LC-MS/MS] 82145 [ELISA] 80101

You are viewing Amphetamines Panel, Blood (Forensic) Test (8215B)

ABCDEFGHIJKLM NOPQRSTUVWXYZ 0123456789 <input type="button" value="SEARCH"/>

SIGN ME UP



Sign me up for e-mail updates