

EXHIBIT 4

LABORATORY NOTEBOOK

No 633

AMGEN, INC.

LIN
EXHIBIT
115

AM 17 009851
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

FRITSCH v. LIN
INTER NOS 102,096
102,097, 102,33A
LIN DOC NO L01097

NOTEBOOK NO. 633
ISSUED TO Dejne
ON 7/10 19 84
DEPARTMENT _____
RETURNED _____ 19 _____

— SCIENTIFIC NOTEBOOK CO. —
5007 WEST DONNA DRIVE
STEVENSVILLE, MICHIGAN 49127

AM 17 009852
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

FRITSCH v. LIN
INTER MDS 102,096
102,097, 102,334
LIN DOC NO L01098

9/19-9/21/84 Western Analysis

Project No. _____

Book No. _____

65

Page No. _____

Large Western analysis - 2 gels to look at the following:

- ① Comparison of B₁₁ Hu EPO B Lot 82 to illustrate similarity in size
ENDO F treatment of 2 as well
 - ② Goldwasser EPO ± ENDO F
 - ③ α Therapeutics EPO ± ENDO F
- ∴ Comparison of size of glycosylated & deglycosyl EPOs from 3 diff. human urinary preps.
- ④ Comparison of ab-reactive EPO in cycles 1 → 9 RB collection + their ENDO F digestion products
 - ⑤ B₁₁ ^{pool for} 18 scale human EPO purification ± ENDO F
 - ⑥ ~~See~~ Cycle ② B₁₁ RB sample ^{is} either treated or not treated ± ENDO F ^{before} ~~then~~ running on gel

⑦ Use α C-terminal peptide in Western ± ENDO F treated EPOs

GEL, TRANSFER + AB PROTOCOL

run 12.5% SDS PAGE - during the day
 - transfer to 0.2 μm Nitrocellulose at 30 v / 0.1 A = 16 hours
 block ± either 10% horse or goat sera in PBS

+bs - ① for part of 1 gel use 208 prep #9 1.05 mg/ml ab at conc of 12 μg/ml - ON then RT in 5% horse sera

Develop ± avidin/biotin reagents - 1.5 hr α mouse Ab, 1.5 hr ABC

② α N₁₆ peptide ab 1:300 diln of 8C196 # 4/25/84 test bleed in 5% goat sera other conditions as in #1

③ α C-terminal peptide ab - 1:200 diln of 8C 278 # 4/25 in 5% goat sera other conditions as in #1

ETOH Precipitation Protocol

from Jeff's Coomassie-stained gel it looked like there was only 1 ENDO F-digestn product + not the usual sublet that we've seen. Difference betw his protocol + ours is that he used ETOH pptn to concentrate sample after ENDO F + before gel loading. The undigested sample was treated the same way.

→ Start ± sample = 15 λ cycle ② + 5 λ ²⁵¹I EPO (23056 cpm). Add 250 μl ^{-80°} ETOH 100%.

→ LE # 12 Sit at -80° for 20 min.

Centrifuge 2' in Eppendorf. Remove sup. See no pellet, but have faith.
 Count pellet + sup (in new tube)

FRITSCH v. LIN
 INTER NOS 102,096
 102,097, 102,334
 LIN DOC NO L01070

To Page No. 66

Read & Understood by me,

Date

Invented by

Date

AM 17 009853

Recorded by

CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

Project No. _____

Book No. _____

TITLE _____

66

From Page No. 65

Sup = 1383 cpm = 6.8% recovery of counts
 pellet = 19061 cpm = 93.2% " " " in pellet.
 20444 total cpm

Let pellet tube stand overnight in hood E fluorescent light on + blowers blowing.
 In AM. all residual EtOH gone

Dissolve in 3x Laemmli sample buffer (LSB) + boil 2' to get into solution.
 20µl

ENDO-F treated sample - cycle (2) EtOH ppt's = SAMPLE # 11
 15µl cycle (2) + 5µl ¹²⁵I-ERD (23,520 cpm) + 20µl ENDO F buffer + 0.5µl ENDO F
 Incubate 5 hrs., 37°
 Cool on ice

Add 500µl 100% EtOH final [EtOH] = 92.93%

Keep at -80° for 20'. Centrifuge + count as above

Sup - 1419 - 6.1% cpm recovered

Pellet - 22004 93.9% cpm recovered.
 23,423

ENDO F. from NEN - Lot #

FRITSCH v. LIN
 INTER NDS 102,096
 102,097, 102,334
 LIN DOC NO 101071

AM 17 009854
 CONFIDENTIAL
 SUBJECT TO PROTECTIVE ORDER

To Page No. 107

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by	

Project No. _____

Book No. _____

67

LE

Page No. 66

REC # 1

for Western #	SAMPLE	Sample Preparation Info.
1	Hi MW BRL prestained MTK = 11 µl	
2	LOT 82 1:20 stock in 0.025% BSA - 6 λ = 52.5 ng	
3	LOT 82 " " - 3 λ = 26.3 ng	
4	CHO Cycle ② - 2 λ = 86.3 ng	
5	CHO Cycle ② - 1 λ = 44.2 ng	
6	LOT 82 + ENDO F - 6 λ ^{LOT} 82 1:20 52.5 ng	① LOT 82 - 10 λ of 1:20 in 0.025% BSA + 40 λ ENDO F BUFFER + 0.5 λ ENDO F = 50.5 λ - incubate 5 hrs, 37° Add 6 µl LSB-10x-boil = 57 µl approx total vol. Use 38 λ for well 6 + 19 λ for well 7
7	LOT 82 + ENDO F - 3 λ ^{LOT} 82 1:20 26.3 ng	
8	CHO Cycle ② - 2 λ + ENDO F 86.3 ng	② CHO cycle ② - 12 λ + 48 λ ENDO F buffer + 1 λ ENDO F = 61 λ TOTAL 5 hr 37° Cool. + 9 λ 10x LSB = 70 λ BOIL TOTAL for 7 λ CHO sample, use 5.5 λ mix. for 2 λ CHO sample, use 11 λ mix
9	CHO Cycle ② - 1 λ + ENDO F 44.2 ng	
10	Goldwasser 1 λ 1:10 for blots = 30 ng	
11	Goldwasser 11 λ 1:10 for blots + ENDO F = 30 ng	③ 33 λ G.C. 1:10 + 66 λ ENDO F buffer - 5 hr / 37° + 0.5 µl ENDO F = 100 λ Add 15 λ 10x LSB = 115 λ - Use 38 λ for wells 18 + 16.
12	dihydrogenated "OLD" - 23 µl = 2.3 U = 33 ng	
13	dihydrogenated "OLD" + ENDO F - 23 µl = 33 ng	④ 23 µl dihydrogenated + 23 µl ENDO F buffer + 0.25 λ ENDO F - 5 hr 37° = 46 µl total vol. + 5 λ 10x LSB - boil - load total sample
14	Lo MW BRL prestained MTK = 11 µl	
15	Goldwasser - 11 λ 1:10 for blots - 30 ng	
16	Goldwasser - 11 λ 1:10 for blots + ENDO F - 30 ng	SEC sample ② above in box
17	CHO - cycle ② - 2 µl = 86.3 ng	
18	CHO - cycle ② - 2 µl + ENDO F 86.3 ng	SEC sample ② above in box
19	CHO cycle ② - 2 µl + ENDO F - 86.3 ng	
20	Lo MW BRL Prestained MTK - 11 µl	

1.6A - ProtA
1.4
1.2 µg/ml

0.0 diln.
1.1 ml db
1.8 & 2.8
1.265

FRITSCH v. LIN
INTER NOS 102,096
102,097, 102,334
LIN DOC NO L01072

AM 17 009855
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

To Page No. 68

Read & Understood by me,

Date

Invented by

Date

Recorded by

Project No. _____
 Book No. _____ TITLE _____

Page No. 61

Used	#	SAMPLE	SAMPLE PREPARATION
	1	BRL 6 MW prestained markers - 11 µl	
	2	Cycle ① - B11 - 1.5 µl = 60 ng	
	3	Cycle ② - B11 - 1.5 µl = 71 ng	
	4	Cycle ③ - B11 - 3 µl = 73 ng	
	5	Cycle ④ - B11 - 4.5 µl = 66 ng	
	6	B11 pool for 1x scale purif = 2 µl = 53 ng	
	7	" " " " - 4 µl = 106 ng	
	8	Ecoli ETO fr. K2 B ⁻ - 2 µl ^{amt ETO unknown}	Make a 1:10 diln. fr. Kris' 20 µl stock + use 2x 2x 1:10 → good control (mad. Priz. is from EK lab. Proj. is in LSA Fridge)
	9	Cycle ① - 1.5 µl = 60 ng + ENDO F	1.5x cycle ① + 10x ^{endo F} buffer + 0.25x ^{5 hrs} ENDO F. 37°, 12x 10x LSA
	10	Cycle ② - 1.5 µl = 71 ng + ENDO F	1.5x cycle ② + 10x " + 0.25x " " + 2x 10x LSA
	11	Cycle ③ - 3 µl = 73 ng + ENDO F	3x cycle ③ + 10x " + 0.25x " " + 2x 10x LSA
	12	Cycle ④ - 4.5 µl = 66 ng + ENDO F	4.5x cycle ④ + 10x " + 0.25x " " + 2x 10x LSA
	13	B11 pool for 1x scale purif = 53 ng + ENDO F	B11 pool 12 µl + 48x ENDO buffer + 0.72 ENDO F = 5 hrs, 37° + 7x 10x LSA = 68x
	14	B11 pool for 1x scale purif = 106 ng + ENDO F	4x 11.3x for the 2x sample + 22.6x for the 4x sample
	15	BRL 6 MW prestained markers - 11 µl	
	16	Cycle ② + ETOH - 2 µl equit = 86 ng	② See description on p. 65-66 Use 2x for the 1x sample Use 4x " " 2x sample
	17	Cycle ③ + ETOH - 2 µl equit = 44 ng	
	18	Cycle ④ + ENDO F + ETOH - 2 µl equit = 86 ng	① See description on p. 65-66 Use 2x for the 1x sample Use 4x for the 2x sample
	19	Cycle ② + ENDO F + ETOH 1 µl equit = 44 ng	
	20	BRL 6 MW prestained markers - 11 µl	

300 dilui.
 196 #, 465 1B
 mit diamino
 unid als

FRITSCH v. LIN
 INTER NOS 102,096
 102,097, 102,334
 LIN DOC NO L01073

AM 17 009856
 CONFIDENTIAL
 SUBJECT TO PROTECTIVE ORDER

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by	


Project No. _____
Book No. _____

TITLE _____

65

From Page No. 68

Get # 1



Conclusions - C10 + lot 82 same size, although C10 very heterogeneous. On + tNDOF, lot 82 gives 1 band + only 1 faint lower MW band, whereas C10 → 3 bands. Intensity lower one is less than higher MW band, but darker than correspond band in lot 82

- Gene's EPD - behaving as usual - glycof = ↓ MW fr. lot 82, although tNDOF digest'n. prod = M_n so faint lower band + 1 def ad. prod.
- dtherapeutics - w same size as C10 + lot 82. Digest'n. prod = lot 82 + lower band is lighter than in C10, but darker than lot 82
- dC terminal sera - light bands - need to repeat - See only top band of digest'n prod. but must repeat E ~ 10x sample load.


To Page No. 70

Witnessed & Understood by me,	Date	Invented by	AM 17 009857	Date	FRITSCH v. LIN
		Recorded by	CONFIDENTIAL		INTER NOS 102,096
			SUBJECT TO PROTECTIVE ORDER		102,097, 102,334
					LIN DOC NO 101024

Project No. _____
 Book No. _____ TITLE _____

Page No. 69

Develop 1:300 Rabbit α NHz Terminus sera





- conclusions
- ① Think there is a tendency towards ↓ MW in ↑ cycle # of ctd samples
 - ② See 3 ENDO F major bands - 2 = weighted doublet
 (why seeing more bands now? 1 lower MW band probably we needed double ab procedure for ENDO F of recombinant material) | degradatiⁿ prod.
 - ③ ETOH pptn. → aggregatiⁿ. + no ↓ in # ENDO F bands
 - ④ E. coli EPD has MW (or migration) even lower than the smallest of the ENDO F fragments of ctd EPD. Why such a discrepancy?

FRITSCH v. LIN 102,097, 102,334
 INTER NOS 102,096 LIN DOC NO L01075

AM 17 009858
 CONFIDENTIAL
 SUBJECT TO PROTECTIVE ORDER



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
	Recorded by		



AM 17 009859
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

FRITSCH v. LIN
INTER NOS 102,091
102,097, 102,334
LIN DDC NO L0107



FRIISCH v. LIN
INTER NOS 102,096
102,097, 102,334
LIN DDC MD L01077

AM 17 009860
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

Project No. _____

Book No. _____

71

TITLE _____

From Page No. _____

[Lined area for transcription or notes]

FRITSCH v. LIN
INTER NOS 102,096
102,097, 102,334
LIN DOC NO L01078

AM 17 009861
CONFIDENTIAL
SUBJECT TO PROTECTIVE ORDER

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by	